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ELITE VEST CAST

Composition

Phosphate bonded, carbon free special investment for model casting - frameworks. Compatible with base- and precious metal alloys (except for titanium). Suitable for speed and overnight pre-heating technique. Indicated for silicone or gel duplication.

Technical data	
- Powder/liquid ratio (Elite Vest Liquid)	100 g / 18-20 ml
- Manual pre-mixing time	15-20"
- Vacuum mixing time (320 rpm)	60"
- Vacuum maintaining time after mixing	10-15"
- Processing time (20-23 °C)	5'
- Time for placing it in the furnace from mixing start	23-25'
- Ideal usage temperature	20-23 °C
- Liquid storage temperature (Elite Vest Liquid)	> 5 °C
- Maximum pre-heating temperature	1050 °C

The indicated times may vary according to material and room temperature. High temperature reduces processing time and makes the material harden faster.

Indications

Suitable for ring and ring-less technique and removable partial dentures made of base- and precious-metal alloys (except titanium).

Classification: UNI EN ISO 15912
Type 2, Class 1,2
Mixing liquid
Elite Vest Liquid 1000ml (universal for Elite Vest investment)

Storage

Store the product in a ventilated place away from light and at suggested temperature of:

- Elite Vest Cast 18-22°C (max 40°C)
- Elite Vest liquid 18-22°C (min 5°C)

Warning: Elite Vest liquid may freeze if exposed to temperatures below 2 °C/36 °F. Once if freezes it can no longer be used.

INSTRUCTIONS FOR USE

Using Elite Vest

Dosage

100 g Elite Vest Cast: 18-20 ml Elite Vest liquid
400 g Elite Vest Cast (bag) 80ml Elite Vest liquid/ silicone duplication – 72ml Elite Vest liquid/ gelatine duplication
We recommend using a graduated measuring cup to dose Elite Vest liquid.

Note: for slow pre-heating, reduce the concentration of the liquid from 5 to 10% according to the alloy and the application.

Mixing ratio

Recommended concentration of Elite Vest liquid (consider the data as a guideline).

% liquid concentration	Application	200 g powder		400 g powder	
		Liquid ml	Distilled water ml	Liquid ml	Distilled water ml
100%	partial dentures, full dentures in one-piece	40	–	80	–
95%		38	2	76	4
90%		36	4	72	8
85%		34	6	68	12
80%	partial dentures, full dentures with attachment	32	8	64	16
75%		30	10	60	20
70%		28	12	56	24
65%		26	14	52	28
60%	Standard partial dentures	24	16	48	32
55%		22	18	44	36
50%		20	20	40	40

% liquid concentration	Application	200 g powder		400 g powder	
		Liquid ml	Distilled water ml	Liquid ml	Distilled water ml
100%	partial dentures, full dentures in one-piece	36	–	72	–
95%		34	2	68	4
90%		32	4	65	7
85%		31	5	61	11
80%	partial dentures, full dentures with attachment	29	7	58	14
75%		27	9	54	18
70%		25	11	50	22
65%		23	13	47	25
60%	Standard partial dentures	22	14	43	29
55%		20	16	40	32
50%		18	18	36	36

Mixing

Dose Elite Vest Cast and Elite Vest Liquid and pour the liquid and then the powder in a bowl; mix manually with a spatula for approximately 15-20". Then stir the mixture under vacuum for 1' using an automatic mixer (we recommend Zhermack Elite Mix). Once the mixture is complete, maintain the vacuum for other 10-15" before using the product. Note: ensure that the mixer creates the vacuum correctly.

Processing

While creating the model investment, we recommend applying Zhermack Tensilab surfactant to the silicone mould before pouring the Elite Vest Cast mixture into it. The minimum time to remove the model investment from the silicone/gel mould is approximately 23-25' at 23°C. Pour the Elite Vest Cast mixture slowly to prevent bubbles. Use a low frequency vibrator to improve the smoothness of the investment.

Speed Pre-heating

Once the mixture is set (23-25' after starting to mix and depending on material and room temperature) roughen the upper surface of the investment before pre-heating . Then place the it in the pre-heated furnace at the final temperature (900-930 °C). Note: gold alloys (Au) with pre-heating temperature ranging between 700 and 750 °C allow for better surface quality results by heating the cylinder investment at 850 - 900 °C first and then reducing the temperature of the furnace down to 700-750 °C.

Time at final temperature: 1 hour

The indicated times refer to only one cylinder. In the event that several cylinders are pre-heated at the same time, we recommend increasing the time by 10' for every additional cylinder.

Overnight Pre-heating

For slow pre-heating, reduce the concentration of the liquid from 5 to 10% according to the alloy and application. Once the mixture is set (23-25' after starting to mix) roughen the upper surface of the investment before pre-heating. Then place it in the furnace at room temperature. Start the multi-stage program with the following parameters

Position	Temperature	Holding time	Heat-rate
1	280 °C	45-60'	3-7 °C (min)
2	580 °C	45-60'	
3	Final temperature	30-45'	

The indicated times depend on the type and number of cylinders in the furnace. Comply with the final temperature recommended for the alloy used.

Maximum pre-heating temperature of the ring: 1050 °C.

IMPORTANT REMARKS:
Elite Vest Plus contains Quartz and Cristobalite; avoid inhaling the powder during use. We recommend wearing a compliant protection mask (FPP 2 – EN149:2001) to prevent damaging the lungs (silicosis or cancer).

GHS08

Remove the powder from the environment only when it is wet. To reduce the amount of powder while deflasking, saturate the cylinder once cold with water. Carry out the sandblasting process in a place equipped with an efficient suction system suitable for fine dust. Dilute Elite Vest liquid using always distilled water. Elite Vest Liquid may freeze if exposed to temperatures below 2 °C/36 °F. Once if freezes it can no longer be used. Avoid procuring the product during the extremely cold season. Clean the investment mixing bowl thoroughly; if they are used also for plaster-based products, they may delay the setting time. Keep the bowls immersed in water between one mixing process and the other. Please use a dedicate mixing bowl to mix investments material. Avoid opening the pre-heating furnace during the wax combustion phase. If you use furnaces with heated chambers at the base, avoid positioning the investment cylinder in direct contact (keep at least 1 cm away).

All the information contained in these instructions for use are indicative and are based on Zhermack's knowledge of the product with respect to the best laboratory results obtained. The combination of Elite Vest with different types of supplementary materials in the various applications may lead to usage differences and different final results.

The warranty covers Elite Vest products for 2 years from the date of manufacture.

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% liquid concentration	Application	200 g powder		400 g powder	
		Liquid ml	Distilled water ml	Liquid ml	Distilled water ml
		40	–	80	–
100%	partial dentures, full dentures in one-piece	38	2	76	4
95%		36	4	72	8
90%		34	6	68	12
85%		32	8	64	16
80%	partial dentures, full dentures with attachment	30	10	60	20
75%		28	12	56	24
65%		26	14	52	28
60%		24	16	48	32
55%	Standard partial dentures	22	18	44	36
50%		20	20	40	40

% liquid concentration	Application	200 g powder		400 g powder	
		Liquid ml	Distilled water ml	Liquid ml	Distilled water ml
		36	–	72	–
100%	partial dentures, full dentures in one-piece	34	2	68	4
95%		32	4	65	7
90%		31	5	61	11
85%		29	7	58	14
80%	partial dentures, full dentures with attachment	27	9	54	18
75%		25	11	50	22
65%		23	13	47	25
60%		22	14	43	29
55%	Standard partial dentures	20	16	40	32
50%		18	18	36	36

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