colorbite rock

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SOLO PER USO DENTALE

Colorbite Rock is a dental vinyl-polysiloxane for bite registration.

Colorbite Rock stands out for its innovative setting time control system, which indicates setting by colour change. Colorbite Rock is lilac during working, and becomes pink once set. The colour changes with the temperature; this allows the clinician to use it safely and easily, according to the correct application times. Colorbite Rock has a pleasant raspberry flavour.

PROPERTIES

Thixotropic. It does not drip.

Very hard after setting (95 Shore A)

Extremely accurate

Unperceivable by the patient, who bites naturally

Dimensionally stable in time and under disinfection

Fast. Time in mouth 45 sec.

Can be repositioned for a check after setting

Easy finishing. It can be finished after setting with a cutting blade or bur.

APPLICATIONS

Bite registrations with all techniques

Intermaxillary centric, protrusion or lateral registration keys

Registration keys for gnathological evaluations

Registrations for diagnostic evaluations of craniomandibular disorders

Orthodontic bite registrations

INSTRUCTIONS FOR USE

Device preparation:

Assemble the device as shown in the drawings (1-2).


Warning: Before assembling the tip, ensure that the two components (base and catalyst) flow out evenly by applying a light pressure on the dispenser lever and extruding a small amount of material, which must be removed (3).

Subsequently, insert the mixing tip and, if required, the intraoral tip in the cartridge (4).

Apply Colorbite Rock according to the desired (occlusal or buccal) technique. If the occlusal technique is used, after extruding the material directly on the teeth, ask the patient to bite. Colorbite Rock is thixotropic; it does not drip, therefore facilitates control. Colorbite Rock is lilac on application; during setting, its colour changes to pink. After setting, the colour must be evenly pink; at that point, the registration can be removed from the oral cavity with no risk of permanent deformation.

The colour change from lilac to pink takes place due to the shift from room temperature (23°C) to oral cavity temperature (35°C). Subsequently, after removal, Colorbite Rock at room temperature reacquires its lilac colour: this change does not in any way alter its properties.

TIMES FOR CLINICAL USE

WARNING: In sensitive subjects, vinylpolysiloxanes may cause irritation or other allergic reactions.

CLEANING AND DISINFECTING THE IMPRESSION

After rinsing well under running water, the impression can be immediately disinfected. Ideal disinfecting is achieved by immersing the impression in **Zeta 7 solution** or spraying it with **Zeta 7 spray** (Zhermack). When using other disinfectants, always follow the manufacturer's instructions.

CLEANING AND DISINFECTING THE DISPENSER

To clean, disinfect or sterilize the dispenser, always follow the manufacturer's instructions..

STORAGE

The registration can be sent to the laboratory like a normal impression, without special precautions, provided that it is protected from mechanical stresses. Store at room temperature.

Colorbite Rock is guaranteed for 36 months if properly stored at temperatures from 5 to 27°C.

WARNINGS

- Hardening of the vinyl polysiloxanes is inhibited by contact with latex gloves. Avoid direct or indirect contact with latex gloves. Wash hands well to eliminate all traces of contaminant or use vinyl gloves.

- Vinyl polysiloxanes are chemically resistant; therefore avoid staining gowns and clothes.

TECHNICAL DATA

Colour	lilac at 23°C (73°F) / pink above 30°C (86°F)
Base/catalyst mixing ratio	1:1
Working time	30 sec.
Time in mouth	45 sec.
Detail reproduction	20 µm
Linear dimensional variation after 24h	- 0.05%
Hardness	95 Shore A

* The times mentioned must be intended from the start of mixing at 23°C . 73°F. Higher temperatures reduce the times, lower temperatures increase them.

IMPORTANT REMARKS: Consulting services on the use of our products, provided in any form, be that verbal, written or through demonstrations, are based upon the current state of dentistry and of our know-how. They are to be considered as non-binding information, also in relation to possible third party claims, and do not exonerate the user from personally checking whether the product is suitable for the planned use. Use and application by the user take place under no possible control by the Company, and are therefore under the user's responsibility. Possible liability for damages is limited to the value of goods supplied by the Company and utilized by the user.

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Zhermack S.p.A. - Italy - via Bovazecchino, 100 - 45021 Badia Polesine (Rovigo) Italy