dispenser D2

1:1 dental material dispenser

IT - EN - FR - ES - DE - PT - NL - FI - EL - DA
NO - SV - SL - RU - PL - CS - TR - RO - HR - HU
SK - LV - LT - ET - BG - ZH - JA - KK

Zhermack S.p.A.

Via Bovazecchino, 100 | 45021 Badia Polesine (RO) Italy

T +39 0425 597611 | F +39 0425 53596

info@zhermack.com | www.zhermack.com

DISASSEMBLING THE DEVICE

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1:1 dental material dispenser.

2. PRODUCT DESCRIPTION

The Dispenser D2 is a non-sterile manual device designed to deliver two-component dental materials (e.g. elastomeric impression materials, dental cements) from a 50 ml cartridge directly inside the patient's oral cavity or into another device (e.g. an impression tray).

The Dispenser D2 is reusable.

3. INDICATIONS FOR USE

This device is intended for use by dental sector professionals and for the purposes of dental treatment. The Dispenser D2 is ideal for delivering two-component dental materials packaged in 50 ml cartridges, such as elastomeric impression materials and dental cements.

4. CONTRAINDICATIONS

There are no known contraindications.

5. SIDE EFFECTS

There are no known side effects.

GENERAL WARNINGS/PRECAUTIONS:

- The product's instructions for use must be kept for the duration of its use.
- This device should only be used by dental professionals within a professional dental environment (public or private health facilities with the necessary legal authorisations).
- To reduce the risk of cross-contamination, always use new, uncontaminated gloves when handling the dispenser.
- **The dispenser is NOT supplied sterile. Before using it for the first time, it must be DISASSEMBLED, CLEANED, DISINFECTED AND STERILISED following the instructions given in sections 6.2.1, 6.2.2, 6.2.3 and 6.2.4.**
- Replace the dispenser if it becomes visibly worn or damaged. The end of the reusable dispenser's life cycle is normally determined by wear and/or damage resulting from use.

6. STEP-BY-STEP INSTRUCTIONS

6.1 PRELIMINARY OPERATIONS

1. Read the instructions prior to use.
2. Disinfect your hands and put on a new pair of disposable gloves.

6.2 CLEANING, DISINFECTING AND STERILISING THE DEVICE

SPECIFIC WARNINGS/PRECAUTIONS:

- The cartridge and the mixing tip must always be removed from the dispenser before the dispenser is cleaned, disinfected and sterilised.

- Also remove the rack from the dispenser prior to cleaning, disinfecting and sterilising the individual components.
- Manual cleaning alone is not sufficient to recondition the dispenser correctly. Cleaning must always be followed by disinfection and sterilisation.
- Manual disinfection alone is not sufficient to recondition the dispenser correctly. Disinfection must always be followed by sterilisation.
- Manually clean and disinfect the dispenser within 1 hour of use.
- Manual disinfection of the dispenser must always be performed using a solution of isopropyl alcohol (<20%) and quaternary ammonium salts (<0.3%) (e.g. CaviWipes®). Follow the chosen disinfectant manufacturer's instructions for contact time.
- The use of disinfectant solutions containing phenolic, iodophoric or organic solutions (e.g. alcohol) may cause stains to appear on the surfaces of the device over time.
- Intermediate level disinfection (inactive against spores) has been validated for dispenser disinfection purposes.
- The dispenser **MUST NOT BE**:
 - sterilised with chemical vapour, dry heat or with cold chemical immersion sterilisers;
 - immersed in disinfectant liquids, as these could cause discolouring or deterioration;
 - cleaned/disinfected in instrument washing machines or thermal disinfection devices;
 - cleaned in ultrasonic baths.
- No automatic methods of disinfection have been validated for the dispenser.
- The dispenser must always be cleaned, disinfected and sterilised before being put away in a place suitable for its storage.
- The dispenser must always be sterilised before use.
- Sterilisation temperature for the dispenser and rack must not exceed 132°C.
- Steam sterilisation in an autoclave has been validated as a method for sterilising the dispenser.
- Once sterilised, the dispenser and rack must be stored in their sterilisation packaging in a closed, clean and dry place well away from splashes, sprays or aerosols of bodily fluids.
- Inspect the sterilisation packaging before reusing the dispenser. If the packaging appears damaged or damp, recondition the dispenser and rack again before use.

6.2.1 DISASSEMBLING THE DEVICE:

a. BEFORE FIRST USE

1. Use a Phillips screwdriver to remove the screws from the two disks securing the rack (fig. 1.1).
2. Remove both disks then press the black lever (fig. 1.2/1) and pull out the rack in the direction shown by the arrow (fig. 1.2/2).
3. Replace the disks on the rack using the screws removed previously (fig. 1.3).
4. The dispenser is now disassembled (fig. 1.4) and ready for cleaning, disinfection and sterilisation (see sections 6.2.2, 6.2.3 and 6.2.4).

b. AFTER EVERY USE

1. Place the contaminated dispenser in a clean, closed container and take it to the area designated for cleaning and disinfection.
2. Remove the disposable cover / barrier envelope.
3. Remove your contaminated gloves, disinfect your hands and put on a pair of clean gloves.
4. Press the black lever on the dispenser (fig. 2.1/1) and pull the rack back to the start of its stroke (fig. 2.1/2).
5. Lift up the cartridge lock (fig. 2.2/1). Remove the cartridge from the dispenser (fig. 2.2/2), then clean and disinfect it following the instructions provided by the cartridge manufacturer.
6. Use a Phillips screwdriver to remove the screws from the two disks securing the rack (fig. 2.3).
7. Remove both disks then press the black lever (fig. 2.4/1) and pull out the rack in the direction shown by the arrow (fig. 2.4/2).
8. Replace the disks on the rack using the screws removed previously (fig. 2.5).
9. The dispenser is now disassembled (fig. 2.6) and ready for cleaning, disinfection and sterilisation (see sections 6.2.2, 6.2.3 and 6.2.4).

6.2.2 MANUAL CLEANING

1. CLEANING THE SURFACES OF THE RACK (component B in figs. 1.4 and 2.6):

- 1.a Using a nonwoven fabric wipe impregnated in a solution of isopropyl alcohol (<20%) and quaternary ammonium salts (<0.3%) (e.g. CaviWipes®), wipe over **all the surfaces of the rack for 30 seconds**. Use additional wipes if dirt remains visible.
- 1.b Brush all the surfaces of the rack for **30 seconds** using a soft bristle brush (e.g. a toothbrush) to eliminate all remaining traces of dirt. If necessary, continue brushing until no dirt remains visible.
- 1.c Wipe all the surfaces of the rack at least **twice (2 times)** to remove all detergent residues, using a clean, non-frayed wipe, dampened (wet but not dripping) with purified water. Use additional clean, damp wipes if necessary.

2. CLEANING THE RECESSES IN THE RACK (component A in figs. 1.4 and 2.6):

- 2.a Using a nonwoven fabric wipe impregnated in a solution of isopropyl alcohol (<20%) and quaternary ammonium salts (<0.3%) (e.g. CaviWipes®), wipe over **all the recesses in the rack for 30 seconds**, passing the wipe to and fro through them. Use additional wipes if dirt remains visible.
- 2.b Brush all the recesses in the rack for **30 seconds** using a soft bristle brush (e.g. a toothbrush) to eliminate all remaining traces of dirt. If necessary, continue brushing until no dirt remains visible.
- 2.c Wipe all the recesses in the rack at least **twice (2 times)** to remove all detergent residues, using a clean, non-frayed wipe, dampened (wet but not dripping) with purified water. Use additional clean, damp wipes if necessary.

3. CLEANING THE SURFACES OF THE DISPENSER (component B in figs. 1.4 and 2.6):

- 3.a Using a nonwoven fabric wipe impregnated in a solution of isopropyl alcohol (<20%) and quaternary ammonium salts (<0.3%) (e.g. CaviWipes®), wipe over **all the surfaces of the dispenser for 30 seconds**. Use additional wipes if dirt remains visible.
- 3.b Brush all the surfaces of the dispenser for **2 minutes and 30 seconds** using a soft bristle brush (e.g. a toothbrush) to eliminate all remaining traces of dirt. If necessary, continue brushing until no dirt remains visible.
- 3.c Wipe all the surfaces of the dispenser at least **twice (2 times)** to remove all detergent residues, using a clean, non-frayed wipe, dampened (wet but not dripping) with purified water. Use additional clean, damp wipes if necessary.

4. CLEANING THE RECESSES IN THE DISPENSER (component B in figs. 1.4 and 2.6):

- 4.a Using a nonwoven fabric wipe impregnated in a solution of isopropyl alcohol (<20%) and quaternary ammonium salts (<0.3%) (e.g. CaviWipes®), wipe over **all the recesses in the dispenser, passing the wipe to and fro through them and into the part that holds the rack for 30 seconds**. Use additional wipes if dirt remains visible.
- 4.b Brush all the recesses in the rack for **30 seconds** using a soft bristle brush (e.g. a toothbrush) to eliminate all remaining traces of dirt. If necessary, continue brushing until no dirt remains visible.
- 4.c Wipe all the recesses in the rack at least **twice (2 times)** to remove all detergent residues, using a clean, non-frayed wipe, dampened (wet but not dripping) with purified water. Use additional clean, damp wipes if necessary.

6.2.3 MANUAL DISINFECTION

1. Using a nonwoven fabric wipe impregnated in a solution of isopropyl alcohol (<20%) and quaternary ammonium salts (<0.3%) (e.g. CaviWipes®), wipe over all the surfaces and recesses of the dispenser and rack, making sure that they **remain visibly wet with the solution for the times specified by the manufacturer of the chosen disinfectant**. If necessary, use additional impregnated wipes to keep all the surfaces and recesses wet for the times specified by the manufacturer of the chosen disinfectant.
2. Wipe all the surfaces and recesses of the dispenser and rack at least **twice (2 times)** to remove all traces of detergent, using a clean, non-frayed wipe, dampened (wet but not dripping) with purified water. Use additional clean, damp wipes if necessary.
3. Leave the dispenser and rack to dry until they are obviously visibly completely dry.
4. Inspect the dispenser and rack to make sure that all visible contamination has been removed. If contaminants remain visible, repeat the steps in sections 6.2.2 and 6.2.3. Replace the device if it becomes discoloured, damaged, worn or misshapen.
5. No additional maintenance or lubrication is recommended.

6.2.4 STERILISATION

1. Package the dispenser and the rack separately in two bags (paper/plastic sterilisation bags) (fig. 3.1) conforming to ISO 11607 and ANSI/AAMI ST76:2017, then place the bags in an autoclave that can perform a pre-vacuum cycle.
2. Use the sterilisation parameters given in the table below:

Parameters	Values
Type of cycle	Pre-vacuum
Reference temperature	132 °C
Exposure time	3 minutes
Drying time	20 minutes

Once sterilised, store the dispenser and rack in their sterilisation packages in a closed, clean and dry place. Inspect the sterilisation packages before assembling the device again. If the packages appear damaged or damp, recondition the dispenser and rack again before use.

6.2.5 ASSEMBLING THE DEVICE

WARNINGS/PRECAUTIONS:

- **The dispenser must ALWAYS be cleaned, disinfected and sterilised before assembly (see sections 6.2.2, 6.2.3 and 6.2.4). This should be done only a short time before the dispenser is reused.**
- To reduce the risk of cross-contamination, always apply disposable covers / barrier envelopes to the dispenser before using it for any new patient, and remove them after every use. The use of disposable covers / barrier envelopes does not eliminate the need for cleaning, disinfection and sterilisation as instructed.

1. Disinfect your hands and put on a new pair of disposable gloves.
2. Remove the sterilised dispenser and rack from their bags (fig. 4.1).
3. Use a clean and sterilised Phillips screwdriver to remove the screws from the two disks on the rack (fig. 4.2).
4. Press the black lever (fig. 4.3/1) and push in the rack in the direction shown by the arrow (fig. 4.3/2).
5. Replace the two disks and their screws on the rack using the same Phillips screwdriver (fig. 4.4).
6. The dispenser is now assembled and ready to receive the cartridge (fig. 4.5).
7. Push in the black lever under the rack (fig. 4.6/1) and pull the rack out from the dispenser to the start of its stroke (fig. 4.6/2)
8. Lift up the cartridge lock (4.7/1).
9. Insert the cartridge, making sure that it is facing the right way. The cartridge must be inserted with the “V” on its base facing downwards (fig. 4.7/2). The cartridge lock can now be lowered easily (fig. 4.8).
10. With the cartridge in the dispenser, proceed as instructed by the cartridge manufacturer.
11. It is advisable to protect the dispenser and cartridge with a disposable cover/barrier envelope (fig. 4.9), to be removed after each use.

7. DISPOSAL

Always wear gloves when handling the device. If the device is contaminated, dispose of it as special waste characterised by a risk of biological contamination. If it is not contaminated, dispose of it according to applicable local regulations.

8. IMPORTANT REMARKS

Information provided in any way, even during demonstrations, does not invalidate the instructions for use. Operators are required to check that the product is suitable for the application envisaged. The manufacturer cannot be held responsible for damage, including to third parties, deriving from failure to follow instructions or from unsuitability for an application. The manufacturer's liability is, in any case, limited to the value of the products supplied. Report any serious incident involving the medical device to the manufacturer and to the relevant authorities.

This symbol identifies a medical device as defined by EU Regulation 2017/745.

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