ELITE TRANSPARENT

ES - MODO DE USO

DESCRIPCIÓN

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EN - INSTRUCTIONS FOR USE

DESCRIPTION

Clear flowing addition cured silicone (vinylpolysiloxane) with medium viscosity and high final hardness (72 Shore A).

USE

Making transparent masks for the moulding of light-curing resins.

ADVANTAGES

The passage of light through Elite Transparent allows complete in-depth light-curing of the composite resins. If the light-curing materials set in absence of oxygen, the dispersion layer is not created on the surface. This results in an increased mechanical resistance and colour stability, as well as in an easier polishing of the prosthesis. Working times are shorter and moulding precision is ensured by the initial fluidity and the high final hardness.

DEVICE PREPARATION:

Assemble the device as shown in the drawings (1-2).

Warning: Before putting on the mixing tip, ensure that the two components (base and catalyst) flow out evenly by applying a light pressure on the dispenser lever and extruding a small amount of material, which must be removed (3).

Subsequently, insert the mixing tip (4).

WARNINGS

Before proceeding to use Elite Transparent, make sure that the surfaces in contact are clean, dry and free of grease. Addition silicones contain a platinum catalyst that can be easily inhibited, thus preventing the product from setting.

According to the adopted technique, the use of Elite Transparent with laboratory silicones (e.g. Platinum – Zhermack), will need a separator to prevent the adhesion between these two materials.

METHOD OF USE

Extrude Elite Transparent leaving the end of the mixing tip immersed in the material to avoid air bubbles. Apply constant pressure to the dispenser and distribute evenly.

At a temperature of 23°C, the working time is 1'30" and setting takes place after 15'.

The impression can then be removed with care and the subsequent phases can be executed without the risk of warping.

TECHNICAL DATA

Working time* (min./sec.)	1'30"
Setting time* (min./sec.)	15'
Shore A hardness after 1 h	≥ 72
Variation in linear dimension (%)	0,185
Tensile strength (N/mm ²)	≥ 6.5 (943 PSI)
Tear strength (N/mm)	≥ 4

*The times mentioned must be intended from the start of mixing at 23°C - 73°F. Higher temperatures reduce the times, lower temperatures increase them.

IMPORTANT REMARKS: Consulting services on the use of our products, provided in any form, be that verbal, written or through demonstrations, are based upon the current state of dentistry and of our know-how. They are to be considered as non-binding information, also in relation to possible third party claims, and do not exonerate the user from personally checking whether the product is suitable for the planned use. Use and application by the user take place under no possible control by the Company, and are therefore under the user's responsibility. Possible liability for damages is limited to the value of goods supplied by the Company and utilized by the user.

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