PLATINUM



USER'S GUIDE

A-Silicone for masks



Platinum, simplify your work

Platinum is a high precision addition silicone available in 95, 85, 85 TOUCH or 75 CAD versions. All silicones in the line are characterised by excellent dimensional stability, short working times and ease of use.

Zhermack offers a complete system of putty silicones with different hardnesses, delivering high performance in numerous applications, including dental aesthetics.

To satisfy the various needs of daily practice, Platinum 75 CAD silicone is designed for scanning with CAD/CAM systems and is ideal in combination with Platinum 95 for excellent reproduction of details. Recommended to make masks and counter-moulds.

INJECTABLE TECHNIQUE WITH PLATINUM 85 TOUCH FOR TEMPORARY RESTORATIONS

Creation of a reinforced temporary restoration with injectable technique, starting with a model prepared on an impression with natural abutments prepared by the dentist.

Materials used: Platinum 85 TOUCH, Acrytemp, Elite Rock.













Unfinished temporary restoration (following mask removal)



INDIRECT TEMPORARY RESTORATIONS BY USING PLATINUM 85 TOUCH

The creation of a temporary restoration enables the dentist to have a functional aesthetic support in the dental practice, before even preparing the abutments in the patient's mouth.

Materials used: Platinum 85 TOUCH, Acrytemp, Elite Rock.









REMOVABLE COMPLETE DENTURE WITH COLD-CURING RESIN FOR POURING TECHNIQUE

The combined use of silicones and cold-curing resins makes it possible to create removable prostheses of high aesthetic and functional quality, saving significant amounts of time compared to the traditional technique which uses heat-curing resins.

Materials used: Platinum 85 TOUCH, Platinum 95, Villacryl SP, Elite Stone.

















FINISHED PROSTHESIS



 Add a central pouring channel if the palate thickness is extremely thin. In this case, pour the resin in from the central channel



Technical features

Produc	t M	lixing time (min:s)	Working time* (min:s)	Setting time* (min:s)	Detail reproduction (µm)	Elastic recovery	Strain in compression	Linear dimensional change (after 24 h)	Hardness (Shore A after 24 h)	Heat resistance
Platinum CAD	75	0:30	1:00	7:00	20	> 99.5 %	< 1 %	0.05 %	75	200 °C
Platinum TOUCH		0:30	1:00	7:00	20	> 99.5 %	< 1 %	0.05 %	85	200 °C
Platinum	85	0:30	2:00	8:00	20	> 99.5 %	< 1 %	0.05 %	85	200 °C
Platinum	95	0:30	2:00	8:00	20	> 99.5 %	< 1 %	0.05 %	95	200 °C

*The times mentioned above are intended from the start of the mixing phase at 23 $^\circ$ C (73 $^\circ$ F).

Codes

Platinum 75 CAD - Scannable A-Silicone



Code		Packaging		
C400741	1 x 800 g Base tub + 1 x 800 g Catalyst tub			

Platinum 85 TOUCH - High precision A-Silicone

	Code	Packaging
	C400750	1 x 4.3 kg Base tub + 1 x 4.3 kg Catalyst tub
	C400751	1 x 200 g Base tub + 1 x 200 g Catalyst tub

Code	Packaging
C400727	1 x 450 g Base tub + 1 x 450 g Catalyst tub
C400725	1 x 1 kg Base tub + 1 x 1 kg Catalyst tub
C400723	1 x 5 kg Base tub + 1 x 5 kg Catalyst tub

Platinum 85 - High precision A-Silicone

Platinum 95 - High precision A-Silicone

	Code	Packaging
	C400720	1 x 450 g Base tub + 1 x 450 g Catalyst tub
	C400700	1 x 1 kg Base tub + 1 x 1 kg Catalyst tub
	C400710	1 x 5 kg Base tub + 1 x 5 kg Catalyst tub



For more information please visit our website: www.zhermack.com



Fulfilling your needs

Zhermack S.p.A. | Via Bovazecchino, 100 | 45021 Badia Polesine (RO) ITALY T +39 0425 597611 | F +39 0425 53596 | info@zhermack.com | www.zhermack.com

