hydrorise implant

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INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

HYDRORISE IMPLANT HEAVY BODY (5:1 mixing ratio): high viscosity vinylpolysiloxane indicated in the single-stage technique with two materials for impression tray filling. To be combined with Hydrorise Implant light body fluid on preparations.

HYDRORISE IMPLANT MEDIUM BODY (5:1 mixing ratio): medium viscosity vinylpolysiloxane indicated in the single-stage technique with one material for impression tray filling and to be injected with syringe on preparations.

HYDRORISE IMPLANT LIGHT BODY (1:1 mixing ratio): low viscosity vinylpolysiloxane indicated in the single-stage technique with two materials for injection on preparations. To be combined with Hydrorise Implant heavy body for impression tray filling.

INDICATIONS FOR USE

Hydrorise Implant is a silicone for taking dental impressions, indicated for single-stage impression technique with one or two viscosities.

USER PROFILE

Hydrorise Implant devices are intended for professional use by dentists.

PACKAGING

380 ml cartridge (Hydrorise Implant heavy body, Hydrorise Implant medium body)

50 ml cartridge (Hydrorise Implant light body)

COMPOSITION

Polyvinilsiloxane Polymer, Methylhydrogensiloxane, Organoplatinum Complex, Surfactant, Pigments, Mint and Vanilla Flavours, Inorganic Fillers.

CONTRAINDICATIONS

No known contraindications.

WARNINGS

Use on patients sensitive to vinylpolysiloxanes is not recommended. Stop using and consult a doctor if irritation, redness or other signs of hypersensitivity occur during use. The product may be irritating to eyes. In case of contact rinse thoroughly with water. Consult a doctor if eye irritation persists. If ingested, mixed impression material can cause digestive problems. Consult a doctor if swallowed.

PRECAUTIONS

Use the product at room temperature 23°C/73°F (higher temperatures reduce working time, lower temperatures increase it). Store the product at a temperature between 5°C and 27°C. Wear compatible gloves, protective goggles, and a mask during use. The use of non-latex dental dams is recommended on patients. Use standard and/or custom impression trays, provided they are not deformable and rigid, preferably with a retention rim. The adhesion of material to the impression tray is essential to avoid errors during extraction from the mouth. To ensure strong bonding, choose the suitable adhesive for the impression tray and for the type of material used. Do not overfill impression trays in order to avoid ingestion and emetic stimulus to patient. Due to the high hardness of the material, it is extremely important to eliminate excessive undercuts before impression taking. Disinfect the impression with a special disinfectant. Prolonged disinfection or disinfection with an unsuitable solution can compromise the impression. After having been rinsed thoroughly under running water, the impression can be immediately disinfected. Follow disinfectant manufacturer's instructions. Do not use product beyond the expiry date on the package. Do not dispose in the environment, but according to local regulations. The product should be stored away from direct sunlight. It is advisable to wrap disposable cartridges and automatic mixers in cellophane after each disinfection stage.

ADVERSE REACTIONS

No adverse reactions have been found.

STEP-BY-STEP INSTRUCTIONS

DEVICE PREPARATION - 380 ml system

- Remove the closing plug by gripping the levers (fig. 5).
- Insert the cartridge in the machine according to the manufacturer's instructions.
- Even out the level of the two components (base and catalyst). Clean, eliminating any excess (fig. 6).
- Insert and correctly fasten the dynamic-static mixing tip (fig. 7).
- Insert the tip lock, fixing correctly the fastening components in the cartridge (fig. 8).
- Start up the machine and wait for material output until its colour becomes uniform and eliminate the non-homogeneous mass (fig. 9).
- The 5:1 Zhermack system cartridge is now ready for use (fig. 10).
- After use, leave the dynamic-static mixing tip inserted.
- To replace it, rotate the tip lock 45° and slide it off (fig. 11). Then remove the dynamic-static mixing tip (fig. 12).
- Follow manufacturer's instructions to remove and replace the cartridge.

WARNINGS

- Always even out material before inserting the dynamic-static mixing tip. Otherwise, start up the machine and wait for material output until its colour becomes uniform, eliminate the non-homogeneous mass and proceed with use.
- Only use original 5:1 Zhermack system cartridges, dynamic-static mixing tips, and tip locks. The use of other mixing tips or syringes for polyether does not guarantee optimal quality and can therefore compromise product characteristics, interfering with the final result.

- If the Hydrorise Implant medium body product is difficult to extrude with the syringe for elastomers, we recommend cutting a few mm off the intraoral tip.

DEVICE PREPARATION - 50 ml system

- Assemble the device as shown in the drawings (fig. 1A-1B-2).
- Before inserting the tip: even out the level of the two components (base : catalyst), pressing down lightly on the dispenser lever, extruding a small amount of material which must be eliminated (fig. 3).
- Insert the mixing tip on the cartridge and the corresponding intraoral tip (fig. 4A-4B).

CLEANSING AND DISINFECTION

After having been rinsed thoroughly under running water, the impression can be immediately disinfected. Follow disinfectant manufacturer's instructions.

SCANNING

Impression scanning can be carried out without the use of opacifying spray using specific scanners. Follow scanner manufacturer's instructions.

RADIOGRAPHY

After the impression has been taken, you can identify any material residues in the oral cavity by means of a normal x-ray check.

MODEL STORAGE AND CASTING

Casting can be performed 90 minutes and up to a maximum of two weeks after disinfection, provided that the impression is kept at room temperature.

DISPENSER CLEANSING AND DISINFECTION

Follow manufacturer's instructions for device cleansing, disinfection and any cold or autoclave sterilization.

WORKING TIME

Hydrorise Implant	Heavy body	Medium body		Light body
	Normal Set	Normal Set	Quick Set	Normal Set
Total working time	2'	2'	1'30"	2'
Time in mouth	3'30"	3'30"	2'00"	3'30"
Total setting time	5'30"	5'30"	3'30"	5'30"

IMPORTANT COMMENTS

Any information provided in any way, including during the demonstrations, does not replace the instructions for use. The operator is responsible for ensuring that the product is suitable for its expected application. The Company shall not be held liable for damages, including those of third parties, which are the consequence of a failure to follow instructions or unsuitability for application and, in any case, within the value limits of the products supplied.

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Cartouche 380 ml (Hydrorise Implant heavy body; Hydrorise Implant medium body)

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