**Zhermack** 

Dental

# acrytemp

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**SOLO PER USO PROFESSIONALE**

### **BIS-ACRYLIC SELF CURING RESIN FOR TEMPORARY CROWNS & BRIDGES**

#### **DESCRIPTION**

Composite two-component resin in a self-mixing cartridge for temporary preparation.

#### **COMPOSITION**

Polyfunctional methacrylic esters.

#### **USES / APPLICATIONS**

Direct or indirect preparation of temporary bridgework and crowns.

Temporary inlay and onlay restoration.

#### **BENEFITS**

Contains acrylates but it is methyl-methacrylate free. Does not exceed a temperature of 40°C (104°F), causing no discomfort to the patient. Optimal resistance to loads and abrasion and excellent dimensional stability.

Aesthetic results are guaranteed by brilliant, stable colours.

Acrytemp is easily repaired in the event of accidental breakage using Acrytemp or a photo-polymerisable flowable composite to restore a perfect chemical bond.

Acrytemp is precise and provides results with exceptionally fine edges.

Its flexibility makes bridgework with long trusses particularly resistant.

#### **INSTRUCTIONS FOR USE**

The cartridge system permits optimum bubble-free mixing, enabling the product to be applied directly on the impression (direct method) or the mask (indirect method).

#### **TAKING AN IMPRESSION**

Before preparation of the stumps, take an impression in silicon or alginate. Use a cutter (such as Putty Cut) to remove areas between teeth and squared off bases to facilitate positioning of the impression in the mouth. If elements are missing in the

back quadrants, make a groove in the impression to create a bar connection between the teeth.

Then proceed with regular clinical procedures (stump preparation, extraction, etc.). If the product is used indirectly in the laboratory, prepare the model as usual.

#### **DISPENSER AND CARTRIDGE**

Pull the rack on the dispenser back to the start of its travel using the black lever located underneath it.

Raise the cartridge block. Put the cartridge in the correct position in the dispenser (the mixing ratio is 4:1).

If it is inserted correctly, the cartridge block will close without difficulty.

Remove the coloured lid. Gently extrude a very small amount of product to obtain precise flow control until the two components flow out together evenly.

Insert the mixer tip, following the V-shaped notch on the coloured ring nut, which must correspond with the appropriate housing on the cartridge. The device is now ready for use.

**NOTE:** the initial extrusion from the mixing cannula (about the size of a pea) should be discarded. Then the following mix will be perfect. This must be done for each new mix.

After each use leave the mixing tip on, it serves as a cap. For more information on filling and using the dispenser, refer to the "Dispenser D2 Instructions" provided.

## APPLICATION

Exert gentle pressure on the dispenser lever to apply Acrytemp in the previously prepared, clean and dry impression.

Position the mixing tip in the deepest part of the impression, leaving it immersed in the material to prevent bubble formation.

The process should take about 50 seconds.

## TEMPORARY PREPARATION

The impression, used as a tray full of Acrytemp, is then repositioned correctly in the mouth or on the model.

After 1-2 minutes at 35°C (95°F) (in the mouth) or 3-4 minutes at a room temperature of 23°C (73.5°F) (on a model) from the start of application, Acrytemp in the elastic phase can be removed without difficulty.

**NOTE:** the setting reaction has to be checked with excess material intra-orally (e.g. with a scaler) or on the model. The oral temperature has a significant effect on the setting reaction and the temporary can only be removed without destruction during the elastic state.

## POST CURING AND FINISHING

Once the temp has been removed from the impression or preparation, eliminate excess material and any squared-off bases. It is best to complete polymerisation by immersing the item in hot water at a temperature of 45°- 55°C (113- 131°F); if this is not possible, wait until 6 minutes have passed since the start of application and then finish using specific rotating tools or very fine abrasives.

**NOTE:** oxygen results in surface inhibition, creating a thin film which can easily be removed with ethyl alcohol or another solvent.

Do not inhale dust produced during finishing. Use adequate protection (a dust suction device or dust mask).

## TEMPORARY CEMENTING

A eugenol-free cement should preferably be used to cement the temp. If a cement containing eugenol is employed, remember that any residues or traces of the material may later give rise to problems involving inhibition of curing reaction of composites or other products, especially acrylic-based products.

## REPAIR

If breakage should occur in spite of the material's excellent mechanical resistance, repair as follows:

- Fracture immediately after temp preparation:  
Fix ends at the point of breakage using more Acrytemp.
- Fracture of an existing, used temp:  
Clean and roughen the surfaces affected and create mechanical retentions.

Then repair the temp using more Acrytemp.

Because such a small quantity of Acrytemp is used in both cases, facilitate polymerisation by immersing the temp in hot water for a few minutes.

**NOTE:** Acrytemp residues which have not yet hardened may be removed using ethyl alcohol or another such solvent.

## PRECAUTIONS

Acrytemp contains acrylates but it is methyl-methacrylate free.

With susceptible patients, sensitization to Acrytemp cannot be excluded. Acrytemp should not be used any more, if allergic reactions are observed. Do not use for patients with allergic reactions against acrylates.

Avoid contact with skin, mucous membrane and eyes.

If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.

### CONSERVATION, STORAGE AND TRANSPORT

The product should be kept at between 2°C (35°F) and 20°C (68°F). Store in a cool dry place (e.g. refrigerator). If kept in the refrigerator, leave the product at room temperature for at least an hour (1h) before use. The product has a shelf life of 24

months from the date of manufacture if kept under this temperature. Do not use Acrytemp beyond its date of expiration.

### CLEANING AND DISINFECTING THE DISPENSER


To clean, disinfect or sterilize the dispenser, always follow the manufacturer's instructions. We recommend **Zeta 3** (Zhermack) disinfectant and **Zeta 2** (Zhermack) sterilisation liquid.

### TECHNICAL DATA

Mixing ratio (Base/Catalyst)	4:1 (automix)
Compressive strength after 24 hours	250 MPa
Flexural strength after 24 hours	65 MPa

### WORKING AND SETTING TIMES

00:00 - 00:50 min	Application in mouth
00:50 - 03:00 min	Setting in mouth (1-2 minutes)
03:00 - 04:30 min	Completion of setting* (at 50 °C)
from 04:30 min	Trimming / polishing

\*If the completion of setting takes place at room temperature, trimming / polishing should not be done before 6 minutes (from start of mixing)

### PRODUCT PRESENTATION

1 50 ml cartridge (automix 4:1), available in colours A1 - A2 - A3 - A3.5 - B1

15 Mixing tips for automix 4:1

### IMPORTANT COMMENTS

Any information provided in any way, including during the demonstrations, does not replace the instructions for use. The operator is responsible for ensuring that the product is suitable for its expected application. The Company shall not be held liable for damages, including those of third parties, which are the consequence of a failure to follow instructions or unsuitability for application and, in any case, within the value limits of the products supplied.

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**ZHERMAK S.p.A.** Via Bovazecchino, 100, 45021 Badia Polesine (Ro), Italy