## Villacryl Repair Free

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**Zhermack**   
Dental

# Villacryl Repair Free

cold-curing acrylic resin  
for denture repairs. Without MMA.

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## Villacryl Repair Free

Cold-curing acrylic resin for denture repairs. Without MMA.

**Product characteristics:** Cold-curing acrylic material for denture repair. Without MMA.

**Classification accordance ISO 20795-1:** Type II, Class I

**Designation:** Repair of complete dentures, partial dentures, acrylic parts in framework dentures and acrylic removable orthodontic appliances.

**Contraindications:** Do not use the product in man sensitized to any of the ingredients.

**Side effects:** Not found in current clinical use of the preparation.

### Color shades:

V2 milk pink veined

V3 dark pink veined

V4 pink veined

T4 pink

### Assortment

Set:

[REF] V320V2Z01 VILLACRYL REPAIR FREE V2 powder 100 g + liquid 50 ml

[REF] V320V3Z01 VILLACRYL REPAIR FREE V3 powder 100 g + liquid 50 ml

[REF] V320V4Z01 VILLACRYL REPAIR FREE V4 powder 100 g + liquid 50 ml

[REF] V320T4Z01 VILLACRYL REPAIR FREE T4 powder 100 g + liquid 50 ml

Powder:

[REF] V320V2P01 VILLACRYL REPAIR FREE V2 powder 100 g

[REF] V320V3P01 VILLACRYL REPAIR FREE V3 powder 100 g

[REF] V320V4P01 VILLACRYL REPAIR FREE V4 powder 100 g

[REF] V320T4P01 VILLACRYL REPAIR FREE T4 powder 100 g

Liquid:

[REF] V320L01 VILLACRYL REPAIR FREE liquid 50 ml

### Chemical composition:

**Powder:** acrylic polymers, benzoyl peroxide, pigments

**Liquid:** ethyl methacrylate, ethylene glycol dimethacrylate, N,N-dimethyl-p-toluidyne

### Method of use:

**Preparation:** Damaged place of the denture have to be cleaned and damped with monomer.

**Recommended following mixing proportions:** 10 g of power and 5 g (5,3 ml) of liquid

**Mixing:** Use a sealable bowl made of chemical-resistant HDPE, ceramic, porcelain or glass to mix the components. Pour the liquid into the bowl first and the powder second. Mix properly with spatula, then leave it in the closed container.

**Laying on:** For hands protection it is recommended to use gloves from polyethylene HDPE.

From time to time check the suitability of the material to be laid on in repair place, especially when the ambient temperature is higher than 20-23°C. In case if dough is fast drying (high ambient temperature) damp the surface of dough placed on the denture (with monomer).

**Polymerization:** Material polymerize in water temperature 50-60°C under pressure of 2 bars for 15 minutes.

**Residual methyl methacrylate monomer:** Material does not contain methyl methacrylate monomer.

**Treatment:** Use standard treatment method. Grind with **Pumice** and polish with **Polishing paste** from Zhermapol.

### Important:

Precisely prepare repaired denture place with standard method.

Keep in the closed vessel mixed powder with liquid,

Time of preparing dough depends on ambient temperature.

The repaired place damp with monomer to ensure proper bonding.

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### Danger and counteraction:

#### Information for dental technicians:

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H335 May cause respiratory irritation.

P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: The place of product contact with skin wash thoroughly with soap and water.

P501 Dispose of contents/container to an approved waste disposal contractor.

While working with the material do not: consume meals and smoke tobacco. Work should be carried out in good ventilated premises. Use stand ventilation. In case of contact of the powder particles with eyes the mechanical irritation may occur. In case of contact acrylic powder with eyes wash it with water and contact physician. Liquid contains ethyl methacrylate - flash point +17,5°C. Avoid direct contact wet acrylic dough with skin. The place of wet acrylic dough contact with skin wash thoroughly with soap water. In case of accidental swallowing contact the doctor. Do not pour liquid into drainage system.

#### Information for dentists:

Possibility to relay infection. Each denture must be cleaned and disinfected before return for repairing.

Denture after repairing procedure must be cleaned and disinfected again before placing in oral cavity.

#### Storage and handling:

Liquid is sensitive to influence of temperature and UV radiation under the influence of which the uncontrolled reaction of polymerization may occur.

Product components shall be stored in sates containers, in dark and good ventilated compartment, with temp. 5-25°C, place not accessible by children.

Do not use the product after expiry date specified on the sales packaging and labels of the producting ingredients.

#### Disposal considerations:

Dispose of contents/container in accordance with applicable laws and regulations. Carton packaging, labels, instructions - the possibility of recycling, dispose of as communal not dangerous waste.

Powder component and a packaging are not dangerous. Dispose of as communal not dangerous waste. In case of plastic container- the possibility of recycling.

Liquid and liquid packaging - hazardous waste. Pass to an authorized waste contractor.

#### NOTE!

#### Fordental use only.

Instructions in written form, oral or information relayed during produce presentation result from actual level of dental technique and technology. It should be treated as obligatory information. Information's obtained from third persons (due to their different level of authorization) do not release the user from personal check whether the preparation is intended to the designated application. Due to the fact that the application takes place without the possibility of control from the producer's side; therefore it limits the producers liability to the value of provided material.

We shall be very grateful for all comments related to products, please mail it to our address: biuro@zhermack.com

**Zhermapol Sp. z o.o.** has established Full Quality Assurance System certified to ISO 9001 / ISO 13485 and under surveillance by **Notified Body Number 0120, SGS United Kingdom Ltd. Systems and Services Certification**. Product has CE mark and comply with Directive 93/42/EEC.

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