
zetaclave B

USER MANUAL
Autoclave

IT | FR | DE | EN | ES

Zhermack 
Dental

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[illegible]

Translation of the original operating instructions**GENERAL INDEX****INTRODUCTION****CHAP. 1**

Page 139	1.1	How to read and use the manual
Page 140	1.2	Terminology and symbols

GENERAL INFORMATIONS**CHAP. 2**

Page 141	2.1	Responsibility
Page 141	2.2	Warranty: general instruction
Page 142	2.3	Authorized personnel
Page 142	2.4	General safety regulations
Page 143	2.5	Information table (annex D - EN13060:2009)
Page 143	2.5.1	Graphic example of various types of cycle

DESCRIPTION OF THE MACHINE**CHAP. 3**

Page 144	3.1	General description of the machine
Page 144	3.2	Table of the components
Page 144	3.2.1	Description of the components
Page 145	3.3	Control panel
Page 145	3.3.1	Keyboard
Page 145	3.3.2	Display
Page 146	3.3.3	Indicators icons
Page 147	3.4	Identification data
Page 147	3.4.1	On-machine labels and symbols

TECHNICAL CHARACTERISTICS**CHAP. 4**

Page 148	4.1	Technical characteristics and operating start
Page 148	4.2	Technical Data
Page 149	4.2.1	Standard accessories
Page 149	4.2.2	Weight and overall dimensions
Page 149	4.2.3	Noise levels

TRANSPORT AND MOUNTING**CHAP. 5**

Page 150	5.1	Warning
Page 150	5.2	Packing and unpacking
Page 150	5.2.1	Unpacking and locating the autoclave
Page 151	5.3	Transport and mounting
Page 151	5.4	Elimination / Dismantling

INSTALLATION AND FIRST-TIME USE**CHAP. 6**

Page 152	6.1	Warnings and precautions
Page 152	6.2	Permitted ambient conditions
Page 152	6.3	Required work space
Page 153	6.3.1	Worktops
Page 153	6.4	Connecting up to the power supply
Page 154	6.5	First installation
Page 154	6.5.1	Water connections for use of machine with demineralizer
Page 155	6.6	Printer
Page 155	6.6.1	Changing the printer roll
Page 156	6.6.2	Print head cleaning
Page 156	6.6.3	Printouts
Page 157	6.7	putting the machine into service
Page 157	6.8	Printer
Page 158	6.9	selecting the language
Page 159	6.10	Purpose of the equipment
Page 159	6.11	Improper use / Forbidden use

USING THE EQUIPMENT**CHAP. 7**

Page 160	7.1	Warning
Page 160	7.2	Safety devices
Page 160	7.3	First time start-up (using the system for the first time)
Page 161	7.3.1	Closure of the door
Page 161	7.3.2	Test cycle
Page 161	7.4	Residual risk

WATER TANK FILLING AND DRAINING**CHAP. 8**

Page 162	8.1	Filling with water
Page 163	8.2	Automatic tank filling with clean water
Page 163	8.3	Manual water tank filling
Page 163	8.4	Filling from water supply mains (Demineralizer) or by means of OSMOSIS system

WATER DISCHARGE**CHAP. 9**

Page 164	9.1	Manual water discharge
Page 165	9.2	Automatic draining of the dirty water tank
Page 165	9.3	Draining and cleaning the clean water tank
Page 165	9.4	Water quality table (DIN EN 285)

STERILISATION TABLE**CHAP. 10**

Page 166	10.1	"B" type autoclave
Page 166	10.2	Night cycle

TEST CYCLES**CHAP. 11**

Page 167	11.1	"BOWIE & DICK" cycle
Page 167	11.2	"HELIX TEST"
Page 167	11.3	VACUUM TEST
Page 168	11.4	Bioassay

RECOMMENDATIONS FOR STERILISATION**CHAP. 12**

Page 168	12.1	Recommendations for sterilisation
----------	-------------	-----------------------------------

ALARMS AND ERRORS**CHAP. 13**

Page 170	13.1	Alarms
Page 172	13.2	Errors

MAINTENANCE**CHAP. 14**

Page 173	14.1	Daily maintenance
Page 173	14.2	Weekly maintenance
Page 173	14.3	Quarterly maintenance
Page 173	14.4	Yearly maintenance
Page 174	14.5	Replacing the door gasket
Page 175	14.6	Unscheduled maintenance
Page 175	14.7	Restarting after a long period of inactivity
Page 175	14.8	Requesting technical assistance

APPENDICES

Page 177	A.1	CE Declaration of conformity
Page 178	A.2	Warranty certificate and identification form
Page 179	A.3	Maintenance and repair worksheet
Page 180	A.4	Authorized service partners
Page 180	A.5	Notes
Page 225		Test certificate

The data given in this manual may be modified without prior notice.

CHAPTER 1: INTRODUCTION

HOW TO READ AND USE THIS MANUAL 1.1

This Use and Maintenance Manual provides information on the utilisation, installation and maintenance of the Zetaclave B medical device.

- Zetaclave B must be used in accordance with the instructions in this manual: it is therefore recommended that you read it carefully prior to installation and start-up. Make sure that you read the entire manual and pay particular attention to messages written in bold type or items highlighted by boxes and/or symbols: the symbols are designed to highlight situations involving danger or requiring caution (see paragraph 1.2).
- Il rispetto delle norme e raccomandazioni riportate nel presente Manuale consente un uso sicuro dell'apparecchiatura ed interventi appropriati sulla stessa.
- This Use and Maintenance Manual is an integral part of the cleaner: it should therefore be kept together with the cleaner so that it can be consulted as and when necessary (keep the booklet in a secure dry place away from sunlight, atmospheric agents, etc.). It must be available for consultation throughout the working life of the medical device even when the machine is sold or finally dismantled.
- It is recommended that you keep this manual constantly updated by integrating any amendments, additions or modifications made by the Manufacturer. Any notes or comments should be made on the blank pages at the rear of the booklet (see Appendices: Notes).
- Use the manual properly so that it is not damaged in any way.
- Do not remove, tear or write on any part of the pages.
- If the manual is lost or damaged and its contents become illegible you can request a replacement copy from the Manufacturer.

The purpose of this Use and Maintenance Manual is to provide users with information on how to operate the machine properly. The manual contains all the information needed to use the device as intended: in particular the manual includes instructions and information on the following:

- correct installation of the machine;
- Detailed description of how the machine works and its components;
- Initial start-up;
- Scheduled maintenance;
- Basic Safety and Accident Prevention info.

This manual contains 14 chapters each of which covers a specific argument.
The appendices at the rear of the manual are a useful complement to these chapters.

We recommend that you read the Instruction Manual carefully before beginning to use the device, to ensure that the operations required are carried out correctly: DO NOT carry out operations other than those described in this booklet. The Manufacturer declines all responsibility for direct or indirect damage to objects, persons or animals deriving from improper use of the appliance.

It is of fundamental importance for the sterilizer and the respective equipment to be used solely for the type of product for which they were designed. We therefore request you **to consult the Declaration of Conformity of this appliance at the end of the booklet: the Class to which the appliance belongs is indicated in the "Category" box.** The "Sterilization Table" (Chap. 10) gives all the information necessary for establishing the type of cycle to be used for sterilizing the various instruments.

1.2 TERMINOLOGY AND SYMBOLS

It is important that you understand the following symbols and their significance as they highlight important information such as critical situations, practical advice, or straightforward information. If you have any doubts as to the significance of a symbol always consult this page.

- Tasks for which failure to observe instructions or tampering with machine parts may put personnel in serious danger are highlighted by the symbol:

DANGER!

This symbol refers to safety standards which must be observed at all times in order to protect both workers and machine.

- Critical situations concerning dangerous areas or tasks in which failure to observe the instructions may render the warranty null and void, or cause damage to machine, connected components or the surrounding area, are highlighted by the symbol:

WARNING!

This symbol refers to safety standards which must be observed carefully in order to guarantee your safety, other people's safety and prevent damage to the machine.

- Any work that fails to comply with instructions in this manual, or any tampering which might automatically render the warranty null and void, is highlighted by the following symbol:

FORBIDDEN!

This symbol refers to actions which must never be carried out (i.e. forbidden actions).

- Any general information and/or advice which may be useful at any time is highlighted by the following symbol:

NOTE!

This symbol highlights information and/or useful advice.

CHAPTER 2: GENERAL INFORMATION

RESPONSABILITY 2.1

Failure to observe the instructions in this Use and Maintenance Manual exonerates the Manufacturer from any liability whatsoever.

For any matters not covered by this manual or about which you have any doubts please contact the Manufacturer directly:

Zhermack S.p.A.
Via Bovazecchino, 100
45021 Badia Polesine • RO • Italy
Tel. +39 0425 597 611 • Fax +39 0425 53 596
<http://www.zhermack.com> e-mail: info@zhermack.com

If machine maintenance work fails to comply with instructions or is done in such a way as to compromise machine integrity or modify its characteristics Zhermack shall be exonerated from any liability as regards worker safety and/or machine performance.

WARRANTY: GENERAL INSTRUCTION 2.2

Included in the Appendices to this Use and Maintenance Manual you will find the Warranty certificate and the ID Form which must be filled out in full.

In general the Warranty is rendered null and void by the following:

- Improper equipment use.
- Incorrect installation.
- Serious negligence of maintenance schedules.
- Modifications to or work on the machine which has not been authorised by the Manufacturer (especially on safety devices).
- Use of non-original spare parts.

IMPORTANT: Always keep the original product packaging throughout the guarantee period. You will need the packaging if the product has to be returned to the Service Centre for repairs or replacement in the event of manufacturing defects or malfunction (with the exception of wilful damage, negligence, failure to perform routine maintenance or damage during transport). If you are not in possession of the original packaging, contact Zhermack and ask for instructions (possible repackaging costs will be charged on delivery). The customer is responsible for all costs associated with unsuitable packaging of the product to be returned to the manufacturer. Zhermack is not liable for faults, malfunctions or damage to the product deriving from or consequent to the use of unsuitable packaging, or packaging not of a quality equal to or better than the packaging used for the original shipment.

Zhermack cannot, as of this moment, be held liable if the machine is modified or tampered with without prior written consent from the Manufacturer. Consequently, any repairs made by unauthorised personnel, use of non-original spare parts or failure to comply with the installation standards specified in this manual shall immediately render the Warranty null and void.

2.3 AUTHORIZED PERSONNEL

Personnel may be divided into the following categories:

• **Operator-user:**

Person trained to use the machines. He/she may carry out all the tasks needed to operate the machines: these tasks include starting/stopping the machine, carrying out routine inspections and any other tasks linked to everyday machine use.

The operator-user must only use the machine when it is fully assembled and the safety devices are working properly as described in this manual.

Before using this machine or carrying out any work on it, the operator must read the entire contents of this manual. This manual has been drawn up to provide the user with essential information on rational, safe utilisation of the machine. The user must comply with the information given at all times.

Furthermore, the user must use the work sheets enclosed with this manual so as to keep a record of maintenance/component replacement work as well as make notes on any machine malfunctions.

• **Zhermack authorised maintenance technician:**

Person authorised to carry out work on the machine under all operating conditions and all safety levels. He/she is also authorised to carry out any mechanical or electrical repairs/adjustments, scheduled maintenance and any component replacement work.

2.4 GENERAL SAFETY REGULATIONS

The operator must read the advice and warnings given below and comply with such regulations at all times: doing so will ensure safe, long-lasting equipment performance:

- The machine must be used by responsible adults only.
- Position the machine in a place where it is not accessible to children.
- Do not use the machine near inflammable or explosive sources.
- Use the machine in dry protected areas.
- Check the condition of the power cable periodically: do not use the appliance if the power cable is not perfectly intact.
- Do not carry out maintenance with the machine running or if it is plugged into the power socket.
- Do not approach the machine with inflammable material.
- Always use personal protection devices, in compliance with the applicable directives.
- Do not use the appliance for purposes other than those mentioned in this Instruction Manual.
- Read the paragraph concerning the technical features carefully before starting up the appliance.

Zhermack cannot be held liable for any damage to persons or objects caused by incorrect maintenance carried out by unqualified personnel or maintenance work that does not comply with the instructions in this manual.

Zhermack is, as of this moment, exonerated from any liability concerning injury to persons or damage to objects that may result from using this equipment in a way other than that described in this manual.

This equipment has electrical parts: therefore, in the event of a fire, no matter how small, use only powder-type extinguishers. Never attempt to put out an electrical fire with water.


INFORMATION TABLE (Annex D - EN13060:2009) 2.5

REQUIREMENTS	B
Dynamic steriliser chamber pressure	X
Air outlet	X
Empty chamber	X
Solid load	X
Small porous objects	X
Small porous loads	X
Full porous load	X
B type hollow load	X
A type hollow load	X
Multiple wrap	X
Drying, solid load	X
Drying, porous load	X
Residual air	

X= available

Graphic example of various types of cycle 2.5.1

type B


N.B. The numbers in the graphs are purely indicative.

CHAPTER 3: DESCRIPTION OF EQUIPMENT

3.1 GENERAL DESCRIPTION OF EQUIPMENT

Zetaclave B is a medical device, devised for steam sterilisation of small instruments and equipment.

3.2 TABLE OF COMPONENTS



3.2.1 Description of the components

1. Control panel	11. Water inlet from mains with use of demineralizer
2. Door microswitch	12. Safety valve
3. Water inlet connector	13. Fuses
4. Tray holder	14. Serial port for STS data logger
5. Bacterial filter	15. Main power cable
6. Dirty water outlet connector	16. Tool for removing tray
7. Printer	17. Emergency funnel
8. Main switch	18. Cleaning sponge
9. Condensate discharge connector	19. Discharge hoses
10. Dirty water discharge with use of demineralizer only	20. Inlet hose with filter

CONTROL PANEL 3.3

In order to understand and use the controls correctly, please refer to paragraphs 3.3.1, 3.3.2 and 3.3.3.

Keyboard 3.3.1**Display 3.3.2**

3.3.3 Indicators icons


*		This icon is displayed shortly after a cycle has started and indicates that the machine door is locked.
		This icon indicates that a cycle is running: you will see the arrows spinning around.
		Icon n° 3 warns that printer is out of paper, or the that printer cover is not properly closed. The machine can execute its cycle without any risk.
		Kind of cycle selected.
		USB icon: the autoclaves are provided with operating software able to communicate with and write to a USB flash drive. Some sterilisers incorporate a USB port, while others can write to a USB flash drive through the special STS Datalogger system connected to the RS232 serial port. The icon is displayed when the USB flash drive is inserted in the connector (or the Datalogger). It is important to read the attached instructions when using either the USB incorporated in the autoclave or the Datalogger system. DATA IS SAVED ONLY AND EXCLUSIVELY IF THE FLASH DRIVE HAS ALREADY BEEN INSERTED IN THE CONNECTOR BEFORE LAUNCHING THE CYCLE.
		This icon warns the user that clean water is at its minimum level: it is not possible to start a new cycle. Fill the machine with water before starting other cycles.
		The clean water tank is full. If the PUMP button is pressed again, the display will show "FULL WATER RESERVOIR", warning the user that it is not possible to continue filling the tank with water.
		While this icon is displayed, no cycle can start: the dirty water tank should be completely drained. Allow all water to drain out before closing the draining tap (see paragraph 9.2).

* see paragraph 3.3.2.

IDENTIFICATION DATA 3.4

The equipment described in this manual has an ID platelet giving information about both the machine and the manufacturer.

- A. Name and address of manufacturer
- B. Equipment code and model
- C. Serial number
- D. Power supply information

The equipment has been built in compliance with the EC directives listed in the enclosed "EC declaration of conformity". When requesting spare parts, information or assistance from your Authorised Servicing Centre always quote the data on the ID plate.

Should the ID plate deteriorate and become illegible, even partially, order a replacement plate from the manufacturer: remember to quote the relevant data.

Do not remove and/or the damage the plate as you must be able to read the data at all times.

**On-equipment warnings 3.4.1**

On the equipment there is a label with the following symbol. It reads:

WARNING! RISK OF ELECTRIC SHOCK. DO NOT OPEN OR REMOVE DOORS OR GUARDS BEFORE DISCONNECTING THE PLUG FROM THE MAINS.

ATTENZIONE! READ CAREFULLY THE INSTRUCTIONS OF THE USER'S MANUAL

WARNING! HIGH HEAT DANGER. YOU MUST USE PROPER TWEEZERS AND THERMAL PROTECTING GLOVES.

A proper earth connection is compulsory. Should it temporarily be necessary to use an extension lead the latter must comply with the standards in force in the country of use.

CHAPTER 4: TECHNICAL CHARACTERISTICS

4.1 TECHNICAL CHARACTERISTICS AND EQUIPMENT OPERATING PRINCIPLE

Zhermack designs and manufactures work instruments that provide cutting-edge technology and first-rate quality. They are specially designed to provide made-to-measure solutions that offer long-lasting performance every single day of their working life.

In making these products Zhermack makes use of the latest design tools, thus ensuring that all finished equipment is as functional as possible. Use of only top-quality materials and thorough testing aimed at providing maximum user safety are an Zhermack constant, making our products safe and internationally competitive.

The autoclave is a device designed for steam sterilization of small sized tools and equipment, and is widely used for medical purposes by general practitioners and dentists, in salons that deal with personal hygiene and body care and also in veterinary surgeries. It is also used for sterilizing materials and equipment that come into contact with blood or physiological liquids, such as instruments used by beauticians, tattoo artists, piercers and hairdressers. The very specific sterilization loads used in these sectors of applications require different performance features for the sterilization cycles.

4.2 TECHNICAL DATA

Mechanical specifications:

Volume	Max 0.14 m ³
Colour	RAL9002
Material	AISI304 Steel/FeP01 Steel

Electrical specifications:

Supply voltage	230 V~ ± 10%
Power	2400 Watt
Frequency	50-60 Hz
Supply cable (L 1.5m)	2+1 x 1,5mm ₂
Fuses	6,3x32mm – T12A
Max. transmitted heat	5,76 MJ/h (1370 Kcal/h)
Insulation class	1

Chamber specifications:

Max. operating pressure	2,3 bar (relative)
Max. operating vacuum	-0,90 bar (relative)
Max. operating temperature	140°C
Material	Stainless steel AISI 304
Dimensions (Ø x P)	18 Lt: 236 x 381,5 mm / 23 Lt: 236 x 530 mm

Clean water tank specifications:

Max. capacity	4,5 litres
Cycles per tankful (before the min. level indicator lights up)	18 Lt= ~6 / 23 Lt= ~5
Material	Polyethylene

Dirty water tank specifications:

Max. capacity	4,5 litres
Cycles per tankful (before the indicator monitoring the max. level of used water lights up)	18 Lt= ~5 / 23 Lt= ~5
Material	Polyethylene
Maximum temperature of drained water	50°C

Bacterial filter data:

Maximum diameter	56 mm
Max. allowed filtering capacity	0.3 micron
No. of cycles before replacement	~300

Tray holder specifications:

Material	Anodised aluminium
----------	--------------------

Dati vassoio:

Material	Anodised aluminium
----------	--------------------

Class of pollution:	2
----------------------------	----------

Standard accessories 4.2.1

The tray holder is supplied together with 4 trays and the following accessories.

Emergency funnel	1	-
Tool for tray extraction and door regulation	1	-
Cleaning sponge	1	-
Fill hose	1	-
Drain tubes	2	-

* Please quote the order code when ordering a spare part. You may order parts from your local Authorised Servicing Centre.

Weight and overall dimensions 4.2.2

DIMENSIONS	18 Lt	23 Lt
LxHxD (type B)	505x400x615 mm	505x400x690 mm
Overall dimensions (door open)	350 mm	
Weight (with empty tanks) (kg)	B= 45 kg	B= 50kg
Weight (with full tanks) (kg)	B= 50 kg	B= 55 kg
Weight per unit area of support	Max 32,46 N/cm ²	

Noise levels 4.2.3

The machine characteristically produces an acoustic pressure level in the operating mode lower than 70 dB (A).

CHAPTER 5: TRANSPORT AND HANDLING

5.1 WARNING


To prevent any injury to persons or damage to things always proceed with the utmost care and attention when mounting the machine. Observe all the precautions and instructions contained in the following paragraphs.

5.2 PACKING AND UNPACKING

The cardboard packaging used for transporting the autoclave IS NOT STERILE.

The sterilizer is packed with all accessories placed inside the chamber. It is wrapped in a protective polyethylene bag and then placed inside a cardboard box. To protect it against accidental impact, it is also padded with polystyrene or cardboard. The box has a wooden frame as further protection during transport.

Keep in a dry and protected area at a temperature of 5-30°C.

IMPORTANT: Always keep the original product packaging throughout the guarantee period. You will need the packaging if the product has to be returned to the Service Centre for repairs or replacement in the event of manufacturing defects or malfunction (with the exception of wilful damage, negligence, failure to perform routine maintenance or damage during transport). If you are not in possession of the original packaging, contact Zhermack and ask for instructions (possible repackaging costs will be charged on delivery). The customer is responsible for all costs associated with unsuitable packaging of the product to be returned to the manufacturer. Zhermack is not liable for faults, malfunctions or damage to the product deriving from or consequent to the use of unsuitable packaging, or packaging not of a quality equal to or better than the packaging used for the original shipment.

5.2.1 Unpacking and locating the autoclave

At least two people are required to remove the equipment from the box, following the instructions given below:

- Open the box and remove the staples to avoid getting scratched or cut while removing the equipment.
- Carefully read the instructions for use.
- Remove the machine from its box by taking hold of the sides, without putting strain on parts made of plastic.
- Place the machine on a perfectly horizontal level surface with a load capacity of at least 70 kg.

- Connect the plug to the mains socket. Do not replace the original plug with other types; do not make additional connections; do not connect to multiple plugs or similar: ensure that the installation to which the sterilizer is connected complies with all relevant legislation and can sustain the specified load (paragraph 4.2).
- Open the door by pushing the DOOR button.
- Take the accessory kit out of the box and switch the equipment off.

Carefully read paragraphs 6.5 and 6.5.1 before starting normal operation.

TRANSPORT AND MOUNTING 5.3

Should it be necessary to move the Zetaclave B proceed as follows:

1. disconnect the equipment from the mains power supply
2. if necessary, empty the water contained in the boiler;
3. hold the equipment by the handle, making sure that you keep it in a vertical position at all times. To this end of, always bear in mind the weight and dimensions of the unit being handled (see section 4.2.2) so as to prevent any damage caused by inobservance of these measurements.
4. pack the Zetaclave B in its original box or, if this is not possible, take all the necessary precautions to protect the steam cleaner against knocks or drops, in that the goods are transported at the owner's risk.

The autoclave is fragile and should therefore be transported with extreme care. DO NOT TURN UPSIDE DOWN.

Failure to observe the above exonerates the Manufacturer from any liability concerning equipment malfunction and consequently renders the Warranty null and void.

When shipping/delivering the equipment to your Authorised Servicing Centre always enclose a copy of the purchase document and a copy of the properly compiled ID form.

If the unit is to be transported by courier, the postal service or rail it is recommended such shipments be insured.

DISPOSAL/DISMANTLING 5.4

When disposing of the **packaging materials** the user must comply with the standards in force in his/her country regarding the following materials:

- Wood/paper: non-polluting materials, but must be recycled properly.
- Polystyrene/plastic: pollutants which must not be burnt (toxic fumes) nor dispersed into the environment, but disposed of in compliance with the standards in force in the country of use.

If the equipment is to be **dismantled** the user must, in compliance with EC directives and the laws in force in the country of use, effect elimination and recycling of the following materials:

- Plastic and glass parts, insulated electrical wiring, rubber parts.
- There are no toxic or corrosive parts.

This product must not be disposed of as household waste but rather, when no longer used, must be collected separately according to UE Directive 2012/19.

CHAPTER 6: INSTALLATION AND FIRST-TIME USE

6.1 WARNINGS AND PRECAUTIONS

Before proceeding with installation make sure that all relevant safety conditions prevail and follow the instructions below carefully.

6.2 PERMITTED AMBIENT CONDITIONS

Unless stated otherwise at the time of order the equipment will be configured to operate properly under the following ambient conditions:

Place of Use	Indoors
Altitude	Up to 2000 m
Working temperature	From 5°C to 40°C
Humidity	Max 80%

Ambient conditions other than those listed above may cause malfunctions or sudden breakdowns. Lighting in the installation area must be sufficient to provide good visibility at every single point on the equipment.

More specifically, luminosity must not be less than 200 lux, lighting must be as uniform as possible and there must be no reflected light as this could dazzle the operator.

The equipment must not be used in workplaces having an atmosphere which is explosive and/or at risk of fire as it has not been designed for use in such areas. Nevertheless, should a fire accidentally break out follow the procedure described in paragraph 2.4.


Work area lighting is an important factor in both personnel safety and the ultimate quality of the work being done. In Italy lighting must comply with a Ministerial Decree Law that clearly defines minimum lighting requirements. In other countries lighting requirements form part of the accident prevention and workplace hygiene standards.

6.3 REQUIRED WORK SPACE

Choosing a good workplace with an appropriate amount of space available for equipment installation is an essential factor in personnel safety, the quality of the work and proper maintenance. This zone must not only be spacious enough to allow for optimum equipment operation but must also be well illuminated, aired, not dusty and not exposed to direct sunlight. Note also that the unit must be positioned so the connection plug can be handled/manipulated with ease.

Worktops 6.3.1

Zetaclave B has been designed for operation on a rigid worktop parallel to the floor. Suitable worktops include service furniture (where the operator is standing) or modelling tables (where work is done sitting down). All worktops must be checked for stability.

CONNECTING UP TO THE POWER SUPPLY 6.4

The user must provide a power connection with relative socket near the machine installation point. The user must also install an adequate electrical circuit breaker upstream from the socket as well as efficacious overload/indirect contact safety devices.

Connection is effected via the safety plug (16A, in compliance with European standards) at the end of the power lead: this must be inserted in the mains socket.

When connecting check:

- that mains voltage and frequency are as indicated on the identification plate (incorrect power supply voltage may damage the machine);
- that the mains supply has a proper, efficient earth connection.

It is forbidden to tamper with the power lead and its plug. When replacing them owing to damage and/or wear please contact an Authorised Service Centre.

A proper earth connection is compulsory. Should it temporarily be necessary to use an extension lead the latter must comply with the standards in force in the country of use.


6.5 FIRST INSTALLATION

Correct installation of the autoclave is key to ensuring its proper operation. Below is an installation checklist:

1. The device must be installed in a laboratory accessible only to authorised personnel.
2. The working environment must be properly lit (see UNI standard 10380) and adequately ventilated.
3. Position the equipment on a perfectly flat and horizontal surface with a minimum load capacity of 70 kg. The autoclave is supplied already levelled. The sterilisation chamber is slightly inclined to the rear of the equipment.
4. Position the autoclave so as to leave sufficient room for inspecting and cleaning the sterilisation chamber.
5. Do not install the autoclave close to washbasins or taps: the cover of this equipment is not watertight.
6. Do not install the equipment close to heat sources (i.e. other autoclaves, ovens or similar).

6.5.1 Water connections for use of machine with demineralizer

The equipment water connections are of primary importance for its correct operation:

**If it is present, remove the silicone hose on "A".
"D" is of nut type, for use with a 4/6mm hose.**

Demineralizer filters must be changed every 3 - 6 months.

To discharge the dirty water with use of the demineralizer ("B"), the tap must always be left open.

NEVER connect "A" to "B" on the same terminal: always use two separate hoses.

If the automatic water filling system ("D") is used, it is COMPULSORY to connect the condensate drainage pipe ("A") to a sewer line and not to the condensate container.

PRINTER 6.6

On starting any cycle, the autoclave starts printing all values concerning the selected cycle, plus the machine model and serial number (paragraph 6.6.3). When the autoclave cycle is completed, the printer stops printing: cut off the printout by using the built-in cutter (simply pull the paper upward).

If the cover is not properly closed, or if the printer is out of paper, the relevant indicator icon will be displayed to warn the user (paragraph 3.3.3).

Refer to the instructions provided by the local health authority for printout-filing.

For a correct and long-lasting storage of the paper printouts, it is necessary to keep them away from light and heat sources.


Changing the printer roll 6.6.1

To change the printer roll: open the cover, insert a new thermal paper roll (maximum width should be 57 mm) in the space provided and lead out the paper over the cover roll while closing the cover. Use only thermal paper. Make sure the thermal paper is correctly oriented: if the paper is inserted the wrong way up, the printouts will be blank.

6.6.2 Print head cleaning

If printouts are difficult to read, the print head should be cleaned using a moistened cloth soaked in alcohol: open the printer cover, remove the paper roll and clean the head (see points shown in the picture). Use compressed air to clear dust from the interior of the printer.


6.6.3 Printouts

How to read and interpret the information provided in the printouts.

Model.	STAR	
Serial number.	RELEASE: 0,02	Software release.
Selected cycle.	S/N: STAR1234	
	DATA: 28:12:04	Cycle start date.
	C/N: 0312	
Info on the cycle underway: minimum cycle temperature and pressure values.	CICLO: BOWIE & DICK	Incremental number of cycles started.
	134° 2,06bar 3,30min	
Column showing time: hour, minutes, seconds.	START	Column showing temperature values in degrees centigrade (°C).
	hh:mm:ss °C bar	Column showing pressure values.
	12:42:25 015,7 0,013	
	12:47:50 025,3 -0,801	
	12:55:41 099,7 0,202	
	12:58:20 068,4 -0,701	
	13:01:45 103,5 0,200	
	13:06:32 068,7 -0,701	
Current cycle phase.	STERILIZZAZIONE	
	13:14:54 135,0 2,207	
	13:15:53 134,9 2,200	
	13:16:52 134,6 2,180	
	13:17:52 134,6 2,179	
	ASCIUGATURA	
	13:18:53 134,2 2,171	
	13:19:51 108,8 0,136	
	13:20:51 095,9 -0,383	
	13:21:50 090,4 -0,126	
Operator's signature.	FINE CICLO: OK	Cycle result.
	OPERATORE:	

PUTTING THE MACHINE INTO SERVICE 6.7

Fig. 1: Place the unit on the horizontal, on the worktop (see section 6.3.1). Connect the plug to the mains: follow the instructions given in section 6.4.

Caution: proceed with care so that the weight of the equipment does not injure your hands. To leave at least 10 cm between the autoclave and the wall.

**SELECTING THE LANGUAGE 6.8**

Fig. 2: Keep SELECT pressed (1). Switch machine on with ON/OFF control (2); Push and hold SELECT until the language appears.


Figure 2

Fig. 3: Press PUMP to change language (1). Switch machine off to confirm choice (2).


Figure 3

6.9 USER MENU

To access the user menu, proceed as described below:

- Switch the machine off.
- Keep SELECT pressed and switch the machine on: release the SELECT button only when the preset language (e.g.: ENGLISH) appears on the display.

Within the user menu, the pushbuttons are used in the following manner:

- START = proceed to the next page.
- DOOR and PUMP = used to modify the settings in the pages.

SELECT = when pressed until the confirmation sound is given, it saves the pre-selected value (not all the pages need to save the pre-selected value).

6.9.1 ENGLISH

use DOOR and PUMP buttons to select the required language. Press the START/STOP button to move on to...

6.9.2 TIME: hh mm ss

This page allows you to set the time. Press the DOOR button to increase the value and the PUMP button to decrease it; use the SELECT button to move the cursor to the value to be changed. Press the START/STOP button to move on to...

6.9.3 DATE: dd mm yy

This page allows you to set the date. Press the DOOR button to increase the value and the PUMP button to decrease it; use the SELECT button to move the cursor to the value to be changed. Press the START/STOP button to move on to...

6.9.5 CHARGE WATER BY...

This page is used to select the water filling mode. If external filling mode is used (demineraliser), set the desired values with the DOOR and PUMP buttons, then hold the SELECT button down for 5 sec. (until the system emits a beep), to save the setting. Press the START/STOP button to move on to...

6.9.6 CYCLES RE-PRINT

This function allows you to save the last cycles executed by the machine with a single action. These cycles can be either printed from the internal or external printer by holding down the SELECT button for 5 sec. Press the START/STOP button to move on to...

6.9.7 EMERGENCY ON/OFF

The "Emergency ON" setting is reserved for specialised technical personnel.

Press the START/STOP button to return to the language setting page.

You can exit the USER MENU at any time by holding down the START/STOP button.

To restore the default settings switch the machine off with the ON/OFF control. On restart press simultaneously SELECT and START/STOP .

PURPOSE OF THE EQUIPMENT 6.10

Zetaclave B is a medical device, devised for steam sterilisation of small instruments and equipment.

IMPROPER USE / FORBIDDEN USE 6.11

Using the equipment for any purpose other than that described in section 6.10 (Purpose of the Equipment) shall be considered improper and thus forbidden.

Improper use renders the Warranty null and void and Zhermack shall not be liable for any damages to objects, workers or third parties.

The main reasons for which the Warranty may be rendered null and void are given in section 2.2 and in the "Warranty Certificate" enclosed with this Manual.

CHAPTER 7: USING THE EQUIPMENT

7.1 WARNINGS

The equipment must not be used with any of its parts disassembled: before using the machine make sure that all its components are properly in place.

For your safety, make sure you pay great attention to the instructions given below.

7.2 SAFETY DEVICES

This equipment includes the following safety devices:

1. Three micro-switches controlling the door and the automatic locking system: they are independent of each other and ensure that the door is closed and properly locked. If problems arise, the user is warned by an alarm and the next cycle cannot start. If any problem arises while the cycle is running, the microprocessor stops the process and immediately reduces the machine pressure level.
2. Two different mechanical thermostats ensure that the temperature of the various components does not accidentally exceed the preset value. Thermostats are manually reset.
3. Four electronic temperature sensors continually monitor all critical points on the machine, to prevent over temperature errors during the operating process.
4. A pressure relief safety valve protects against the risk of explosion.
5. An electronic pressure transducer checks all the solenoid valves and opens them if an overpressure is detected.

7.3 FIRST TIME START-UP (USING THE SYSTEM FOR THE FIRST TIME)

After installing the autoclave (as per paragraphs 6.5 and 6.5.1), make sure that all water connections are watertight. Proceed as detailed below:

1. Start the equipment using the main switch.
2. Enter the USER MENU (paragraph 6.9) and configure the system in accordance with the selected filling mode (if water filling is from the water mains, check that inlet water connections are as per paragraph 6.5.1).

The machine is delivered without any water inside

3. Fill the tanks with water: Once the clean water LIV MAX (MAX. LEVEL) icon is displayed, the test cycle can start (paragraph 7.3.2).
4. Press DOOR to open: insert at least the tray holder and start the cycle by using the START button: the autoclave will run automatically until the end of the cycle, when a warning sound will be heard and the display will show "END CYCLE".

Closure of the door 7.3.1

Keep the door pushed in until the display shows the door locked icon.

To close the door, it is NOT necessary to press the START button as the motor starts automatically. To start the cycle, press the START button once the door is closed.

Test cycle 7.3.2

The purpose of the test cycle is to ensure that the autoclave is in working order, has not suffered any damage during transport and is not subject to any operational faults due to technical causes. It is recommended that this test is based on BOWIE & DICK cycle.

- Close the door and press the button START (please to read point 7.3.1)
 - BOWIE & DICK: While the cycle is running, the display will show the temperature – pressure– time remaining– selected cycle number – type of cycle - current cycle phase - any warning icons. Insert the tray holder and the trays inside the machine as indicated at paragraph 11.1. The vacuum pump starts removing air from the chamber until the preset value is reached (VACUUM 1), after which the rise phase will start (PREHEATING). During this phase a slight buzz coming from the water injection pump (used to inject steam into the chamber) can be heard. Once the preset vacuum degree is reached the EXHAUST phase will start; on completion, another two vacuum phases and two preheating phases will be carried out until the cycle preset values are reached (paragraph 10). The STERILISATION phase will now start. during the exposure period, pressure and temperature will be continually software-controlled to ensure that the sterilisation is successful. Alarms will trip in the event of problems (see paragraphs 13 and 13.1).
- On completion of the sterilisation phase, the DRYING phase will begin: pressure inside the chamber will be released and the final vacuum phase will start (the pump will suck any steam to considerably improve the final drying of the sterilised instruments). The door can only be opened when the display shows END CYCLE. In the event of residual positive or negative pressure, the door lock system will not open (proceed as indicated in paragraph 3.3.1, DOOR and PUMP).
- Remove the instruments by using the special handle supplied with the equipment: use protective gloves to avoid getting burnt.

RESIDUAL RISK 7.4

Residual risk is:

- a danger that cannot be completely eliminated through design and application of safety devices.
- a potential danger which cannot be highlighted.

To prevent accidents always observe the instructions in this manual. Should you have any doubts at all please contact the Manufacturer or your Authorised Servicing Centre.

The above-mentioned dangers are indicated by special labels attached to the equipment, as described in sub-section 3.4.1.

CHAPTER 8: WATER TANK FILLING

8.1 FILLING WITH WATER


Figure 4

Fig. 4: Insert the end of the hose without the filter into the connector. The end with the filter must be inserted into the water jerry-can.

5-litre jerry-can (full load: 4.5 l).



Figure 5

Fig. 5: Press PUMP to fill.
The required filling time is shown on the display.


Figure 6

Fig. 6: Detach the filling hose and insert the tray holders in the chamber.
Press SELECT to select a programme.


Figure 7


Fig. 7: Keep the door closed and press START/STOP.

AUTOMATIC TANK FILLING WITH CLEAN WATER 8.2

The total capacity of the clean water tank is 4.5 litres. In order to fill the tank, connect the applied hose to the hose connection fitting, at top, left-hand side of the chamber and press the PUMP button: the filling pump will run for a maximum of 220 sec. When the icon signalling the maximum level (paragraph 3.3.3) lights up on the display, the water supply will continue for a further 5 sec., after which the pump will stop automatically. The filling operation can be stopped at any time by pressing the PUMP button. If the water is at minimum level, as indicated by the relevant icon on the display, the cycle cannot start.

MANUAL WATER TANK FILLING 8.3

In case of failure of the filling pump, the tank can be filled manually. Start the machine and unscrew the cap located alongside the display, then, using the funnel supplied, fill the tank as shown in the photo: When the display shows the maximum level icon (paragraph 3.3.3), manual water filling must be stopped.

Use only demineralised water (water quality is detailed in paragraph 9.5).

Funnel is supplied as an accessory: remove the cap and fit the silicone hose connected to the funnel.

The filling operation should be carried out with the machine ON: this makes water filling easier to control.

On completion of operation screw the cap back in place by finger tight.

FILLING FROM WATER SUPPLY MAINS (DEMINERALIZER) OR BY MEANS OF OSMOSIS SYSTEM 8.4

If a water purification system (Demineralizer or OSMOSIS system) is present, the water connections as described at point 6.5.1 must be used and the software must be programmed with the pre-selected setting. Proceed as described below:

- Switch the machine off.
- Keep SELECT pressed and switch the machine on: release the SELECT button only when the preset language (e.g.: ENGLISH) appears on the display.
- Press START to reach the page in which it is possible to select the clean water tank filling mode. The various possibilities are:
 - **FILLING BY MEANS OF A PUMP** (the machine takes in the water from the front: point 8.2).
 - **FILLING FROM THE WATER SUPPLY MAINS** (the machine requires a purification system with resins).
 - **OSMOSIS** (the machine needs an osmotic purification system).
- Select the pre-selected method with PUMP, press SELECT until a sound confirms that the selected value is saved in memory, then switch the machine off.
- Switch the machine on and check all the hydraulic connections during the tank filling.

Read the user instructions for the purification system to be installed (Demineraliser or OSMOSIS system) and follow the directions given. The manufacturer declines all responsibility for problems deriving from installation and connections carried out by unauthorized personnel.

After installing the purification system, switch the machine on and check all the hydraulic connections during tank filling.

CHAPTER 9: WATER DISCHARGE

9.1 MANUAL WATER DISCHARGE


Figure 8

Fig. 8: Maximum dirty water level symbol.

The capacity of the dirty water tank is 4.5 litres.


Figure 9

Fig. 9: Press DOOR to open.

Dirty water may contain contaminated residues and it is therefore recommended that protective gloves are used while draining water from the equipment.


Figure 10

Fig. 10: Connect the hose (supplied) to the discharge connector. Push the connector towards the autoclave and discharge all the water.

Wait for all the water to go out before disconnecting the tube.


Figure 11

Fig. 11: Pull the connector and remove the hose once the task has been completed.

Bacterial filter: replace after approx. 300 cycles or every three months at the latest.

AUTOMATIC DRAINING OF THE DIRTY WATER TANK 9.2

To eliminate the need to empty the dirty water tank manually, use the tap at the rear of the autoclave. Connect a hose to the tap, secure it with a metal hose clamp, and insert then the other end of the hose into a container or connect it directly to a water drain (paragraph 6.5.1). Leave the tap open: the dirty water will be drained automatically into the collection container or into the direct drain.

DRAINING AND CLEANING THE CLEAN WATER TANK 9.3

In order to empty the tank and gain access for cleaning purposes, proceed as indicated in paragraph 6.9.8. This operation should be performed only after disconnecting the water supply line from the mains (if the sterilizer is connected to an external water purifier).

WATER QUALITY TABLE (DIN EN 285) 9.4

CEN STANDARD DIN EN 285	Max value
Evaporation residue	10 mg/l
Silicon oxide (SiO ₂)	1 mg/l
Iron	0,2 mg/l
Cadmium	0,005 mg/l
Lead	0,05 mg/l
Heavy metal residues (except iron, cadmium and lead)	0,1 mg/l
Chloride	2 mg/l
Phosphate	0,5 mg/l
Conductivity (at 20°C)	15 µs/cm
pH value	da 5 a 7
Look	colourless, clean, without sediment
Hardness	0,02 mmol/l

The use of water containing contaminants at concentrations exceeding those given in the table can considerably shorten the operating life of the equipment, seriously damaging its components, in particular way to the vaporizer, and invalidating the manufacturer's warranty.

CHAPTER 10: STERILISATION TABLE

Sharp surgical instruments must be wrapped to ensure they are sterile at time of use. Data shown in table are only approximate: the type of sterilisation cycle selected should be based on data supplied by the manufacturer of the items to be sterilised. This autoclave is not suitable for the sterilisation of liquids. The overall cycle time can change due to various factors (e.g. weight and type of load, etc.).

To consult the Declaration of Conformity of this appliance at the end of the booklet: the Class to which the appliance belongs is indicated in the "Category" box.

10.1 "B" TYPE AUTOCLAVE

Materials and instruments to be sterilised	Programme	Vacuum phases	Total time*
Fragile hollow** and stainless hollow (in packets and not in packets)	121° Hollow wrapped	4	26 min
Stainless hollow (in packets)	134° Hollow wrapped	4	14 min
Rubber and fragile solids (in packets)	121° Wrapped solids	2	26 min
Stainless solids (in packets)	134° Wrapped solids	2	14 min
Solids and stainless hollow (in packets)	134° Prion	4	30 min
Delicate, hollow and porous (in packets)	121° Porous	4	31 min
Inoxidizable solid, hollow and porous (in packets)	134° Porous	4	19 min
Rubber and fragile solids (not in packets)	121° Fast	2	20 min
Stainless solids (not in packets)	134° Fast	2	8 min
Hollow in rubber and inoxidizable (not in packets)	134° Hollow open	4	8 min
Test cycle for Helix and Bowie & Dick test	134° Helix - Bowie & Dick test	4	7,5 min
Vacuum test cycle	<40° Vacuum test	1	15 min

*: Drying stage included, preheat excluded. Preheat may vary from 20 to 25 minutes depending on the load. For the Zetaclave B 23l the time is increased by 5 minutes. The tolerable load can increase of around 1 kg (0.5 kg on "Porous" cycles).

**: Hollow instruments of type A and B as per EN13060.

To cancel a programme while running, press and hold the START/STOP button for at least two seconds. When the autoclave displays alarm code AL0001, perform a reset by using the SELECT - START/STOP buttons (as described in chapters 13 and 13.1).

10.2 NIGHT CYCLE

If the autoclave is not in use, it will switch to energy-saving mode with only the display background lit. Pressing any key (except the START/STOP button), the display will show the result of the last operation carried out (e.g. END CYCLE). Any cycle can be made a "Night cycle".

After the "night cycle", to the opening of the door, it is normal to find it condenses of water on the gasket of the porthole and on the fund of the sterilization chamber. In case of alarms repeat the operation!

CHAPTER 11: TEST CYCLES

"BOWIE & DICK" CYCLE 11.1

The Bowie & Dick cycle verifies that the steam penetrates properly inside a porous load. Carry out the test cycle after removing all trays from the autoclave chamber but the central one: place the test pack without any other instrument on this tray. Select the "Helix / B&D" cycle and start the cycle. The cycle result is confirmed by the test-pack result.

"HELIX TEST" CYCLE 11.2

This test should be performed to check steam penetration into a hollow load. Carry out the test cycle after removing all trays from the autoclave chamber but the central one: place the "Helix-Indicator System" on this tray (without any other instrument). Select the "Helix / B&D" cycle and start the cycle. The cycle result is confirmed by the test paper included in the test material (see drawing below).

Bowie & Dick test: positive

Helix test: positive / negative

VACUUM TEST 11.3

The Vacuum Test cycle is used to detect any pressure losses in the sterilisation chamber. This test should be carried out on an empty machine, prior to starting other sterilisation cycles. **The cycle cannot start if the temperature inside the sterilisation temperature is >40°.** Select "Vacuum Test" and start the cycle. The autoclave reaches the selected vacuum degree and will hold it for 15 minutes. The test result is given by the "END CYCLE" message appearing on screen and by the relevant printout. If the machine failed the test (AL0600 or AL0601), the door seal should be checked, cleaned or, if necessary, replaced (paragraph 14.1). Similarly, the sterilisation chamber rim should be checked also. Finally, the test should be repeated. A failed Vacuum Test does not bar the use of the steriliser in the near future. However, it is recommended that you contact a technical support centre as sterilisation cycles may be affected in the long term.

11.4 BIOASSAY

Together with other chemical tests, a bioassay may be required. This assay consists in sterilising one or more vials containing biological spores, together with the normal items to be sterilised. On completion of the cycle, remove the vials and leave them to cool for a few minutes (following the manufacturer's instructions for the control procedure). Sterilised vials should normally be broken using the special tools supplied by the manufacturer and then inserted into a specific incubator: together with these, another vial that has not undergone sterilisation should be added for comparison. After the incubation period, the difference in colour in the sterilised vials will indicate whether the cycle was successful.

CHAPTER 12: RECOMMENDATIONS FOR STERILISATION

12.1 RECOMMENDATIONS FOR STERILISATION

In order to prolong the life of both autoclave components and instruments, recommended procedures should be followed as well as local health authority instructions. Below is a list of some precautions to be followed.

1. Instruments should be cleaned with appropriate disinfectants immediately after use.
2. Brush the instruments to remove any residue.
3. Rinse the instruments in running water at room temperature.
4. Submit the instruments to ultrasonic bath treatment.
5. Rinse the instruments in demineralised water at room temperature.
6. Dry the instruments thoroughly.
7. Place the instruments on the sterilizer tray, so that the bags do not overlap. If you need to sterilise instruments which are not wrapped in bags, the tray should be covered with suitable napkins, to ensure that each sterilised instrument is perfectly dried. Follow the manufacturers' instructions for each instrument.
8. Instruments such as scissors or forceps should be opened slightly.
Mirrors should be positioned face down.
9. Place bags with the paper side up.
10. If empty containers are sterilised, place them upside down to prevent water from accumulating.

The instructions listed above show how important the correct preparation of instruments is to successful sterilisation. Even a single instrument carrying traces of disinfectant placed in the steriliser may cause damage to the sterilisation chamber and the instruments inside. In such a case, the sterilisation process might be adversely affected even if no alarm is triggered.


CHAPTER 13: ALARMS AND ERRORS

If alarms are visible on the display (paragraph 13.1), these will prevent any subsequent operation: **in order to reset an alarm press START and SELECT buttons at the same time** until the display temporarily goes out. Alarms are also logged on the printout (see table below). Conversely, errors (paragraph 13.2) do not allow the cycle to restart but indicate that an operation is required before continuing with sterilisation (e.g.: "DRAIN WATER").

In the event of alarms, the cycle running should be considered invalid (i.e. material not sterilised).

Interpreting ALARM codes:

13.1 ALARMS

CODE AND MEANING	PROBLEM	SOLUTION Reset alarm = START/STOP + SELECT for 5 sec.
AL0001 Cycle stopped by user.	Occurs if the START/STOP button is pressed for 1 second or longer	Reset the alarm and repeat the cycle.
AL0002 No mains supply	Caused by voltage peaks higher than 10% or by power failures.	
AL0003 Door open during cycle	Occurs if a control micro-switch detects an open door while the cycle is running.	Reset the alarm and repeat the cycle. If the problem persists, contact technical support.
AL0004 Timer stopped	The timer will stop when PCB battery power is insufficient.	Reset the time and date again following the instructions in paragraph 6.9.2 and 6.9.3). Leave the autoclave ON for at least one hour.
AL0011 No 1st vacuum	The alarm is triggered only if the first vacuum is not reached.	Reset the alarm, clean the door gasket as described in point 14.1, dry the instruments to be sterilised, make sure the pipes are not obstructed by bags in the chamber and then repeat the cycle.
AL0012 No 2nd vacuum	The alarm is triggered only if the second vacuum is not reached.	
AL0013 No 3rd vacuum	The alarm is triggered only if the third vacuum is not reached.	
AL0015 No vacuum when drying	Occurs if the preset vacuum degree is not reached during drying. If this alarm is triggered, the machine has already completed the sterilisation.	Reset the alarm and repeat the cycle. Check if the problem persists in the following cycles: in this case contact technical support.
AL0021 No 1st rise	The machine is unable to reach the 1st pressure value set.	Reset the alarm, fill the tank until the maximum icon appears on screen and repeat the cycle: if the problem persists in the following cycles, contact technical support.
AL0022 No 2nd rise	The machine is unable to reach the 2nd pressure value set.	
AL0024 No final rise	The machine is unable to reach the operating pressure.	
AL0031 No 1st pressure release	After reaching the first set pressure, the machine triggers an alarm.	Reset the alarm, remove the tray holder, clean the sterilisation chamber inside and repeat the cycle.
AL0032 No 2nd pressure release	After reaching the second set pressure, the machine triggers an alarm.	
AL0034 No final pressure release	The machine fails to release pressure while drying.	
AL0100 T1 probe encoding error	Alarm due to PCB self-diagnosis.	Reset the alarm, switch the machine off and on: if the problem persists, contact technical support.
AL0101 OPEN T1	System detects that T1 probe is open.	
AL0102 T1 short circuit	System detects short-circuit in T1 probe.	
AL0110 High temp. in T1 probe	T1 probe has exceeded the preset cycle temperature.	Reset the alarm then wait for 10 minutes with the door open. Repeat the cycle: if the problem continues contact the after sales service.
AL0111 Low temp. in T1 probe during sterilisation	Temperature detected by T1 probe during sterilisation is below the minimum set value.	
AL0200 T2 probe encoding error	Alarm due to PCB self-diagnosis.	Reset the alarm, switch the machine off and on: if the problem persists, contact technical support.
AL0201 OPEN T2	System detects that T2 probe is open.	
AL0202 T2 short circuit	System detects short-circuit in T2 probe.	

CODE AND MEANING	PROBLEM	SOLUTION Reset alarm = START/STOP + SELECT for 5 sec.
AL0210 High temp. in T2 probe	T2 probe has exceeded the preset cycle temperature.	Reset the alarm then wait for 10 minutes with the door open. Repeat the cycle: if the problem continues contact the after sales service.
AL0211 Low temp. in T2 probe during sterilisation	T2 probe detects excessively low values during sterilisation.	
AL0300 P probe encoding error	Alarm due to PCB self-diagnosis.	Reset the alarm, switch the machine off and on: if the problem persists, contact technical support.
AL0301 OPEN P	System detects that P probe is open.	
AL0302 P short circuit	System detects short-circuit in P probe.	
AL0310 High pressure during sterilisation	P probe detects excessively high values during sterilisation.	Reset the alarm, remove the tray holder, clean the sterilisation chamber inside and repeat the cycle.
AL0311 Low pressure during sterilisation	P probe detects excessively low values during sterilisation.	Reset the alarm, clean the gasket and repeat the cycle: if the problem continues call the after sales service.
AL0400 TVP probe encoding error	Alarm due to PCB self-diagnosis.	Reset the alarm, switch the machine off and on: if the problem persists, contact technical support.
AL0401 OPEN TVP	System detects that TVP probe is open.	
AL0402 TVP short circuit	System detects short-circuit in TVP probe.	
AL0404 Low temp. in TVP probe	TVP probe does not reach operating temperature.	Reset the alarm and repeat the cycle. If the problem persists, contact technical support.
AL0405 High temp. in TVP probe	TVP probe has reached a temperature higher than operating threshold.	Reset the alarm then wait for 10 minutes with the door open. Repeat the cycle: if the problem continues contact the after sales service.
AL0500 TRS probe encoding error	Alarm due to PCB self-diagnosis.	Reset the alarm, switch the machine off and on: if the problem persists, contact technical support.
AL0501 OPEN TRS	System detects that TRS probe is open.	
AL0502 TRS short circuit	System detects short-circuit in TRS probe.	
AL0504 Low temp. in TRS probe	Temperature does not reach operating temperature.	Reset the alarm and repeat the cycle. If the problem persists, contact technical support.
AL0505 High temp. in TRS probe	TRS probe has reached a temperature higher than operating threshold.	Reset the alarm then wait for 10 minutes with the door open. Repeat the cycle: if the problem continues contact the after sales service.
AL0600 Excessive drop during the VACUUM TEST stabilisation phase	An excessive pressure drop occurred during the first 5 minutes of the stabilisation phase, while performing the VACUUM TEST cycle.	Take the tray holder out: clean and dry thoroughly the sterilisation chamber. Repeat the cycle: if the problem continues you have to do a Vacuum test cycle when the autoclave is cold.
AL0601 Excessive drop during the VACUUM TEST maintenance phase	An excessive pressure drop occurred during the first 10 minutes of the maintenance phase, while performing the VACUUM TEST cycle.	
AL0700 T1/T2 comparison	Conflicting temperature values detected by the two interior sensors during sterilisation.	

13.2 ERRORS

The following table contains a list of messages which may appear on display:

MESSAGE	CAUSE	SOLUTION
OPEN THE DOOR	At the time of power-up you are prompted to open the door so that the machine can perform a pressure check-up.	Open the door to allow automatic setting of atmospheric pressure (two beeps).
DOOR OPEN	A cycle has been started with the door open.	Close the door and start the cycle.
DOOR NOT LOCKED	Even though the door is perfectly closed, it is read as "not locked" due to mechanical reasons (if this error persists, contact technical support).	Press the DOOR button, open and close the door, then press the START button to start the cycle.
TO LOAD WATER	A cycle has been started with the clean water minimum level message flashing on the display.	Fill the tank with clean water (paragraph 8.2).
TO DRAIN WATER	A cycle has been started with the dirty water maximum level message flashing on the display.	Drain dirty water from the tank (paragraph 9.2).
FULL WATER RESERVOIR	Attempt to fill the tank with clean water with the full tank icon already present on the display.	Disconnect the water hose. It is now possible to start a new cycle.
DOOR STILL LOCKED	Door lock system fails to open completely when cycle terminates.	Close the door again and start a new cycle: stop the cycle after a few seconds with the START/STOP button (paragraph 3.3.1), reset the alarm (paragraph 13), then try opening the door with the DOOR button. If the error persists, contact technical support.
CHAMBER TEMP.>40°	Attempt to start a VACUUM TEST cycle with a chamber temperature higher than 40°C; wait for temperatures to drop before starting the test cycle.	Wait for the temperature on display to go down below 40° with the door open; then the test cycle can be started
REPLACE RESINS	The autoclave cycle counter prompts the operator to check water purifying system resin quality.	If everything is OK, keep the "LOAD WATER" button pressed for 5" (until the sound sets off).
SERVICE	The machine has run an excessively large number of cycles and needs periodic overhauling by a technician.	Contact your dealer as soon as possible and request the services of a technician.

CHAPTER 14: MAINTENANCE

Before cleaning the equipment exterior always turn the main switch to OFF. It is strictly forbidden to clean the equipment while it is connected to the mains power supply.

DAILY MAINTENANCE 14.1

Daily maintenance includes keeping the door seals in good working order, cleaning the chamber rim area (very important for a successful test cycle) and check the water levels in the tanks.

- **DOOR SEAL:** clean the door seal using the soft part of the sponge supplied with the equipment. Cleaning should carry out to remove any impurities which might affect the test cycles.
- **CHAMBER RIM:** it is the outer rim of the sterilisation chamber with which the seal makes contact. Use the coarse part of the sponge supplied.
- **WATER LEVELS** (paragraph 3.3.3): before starting another sterilisation cycle, check the water level in the tanks.
- **GENERAL SURFACE CLEANING:** use a cloth to remove dust and other deposits from the top of the machine.

WEEKLY MAINTENANCE 14.2

Weekly maintenance requires a visual check and cleaning of the inside of the sterilisation chamber. Remove the trays and tray holder from the chamber before cleaning.

- **CHAMBER INTERIOR:** Use the coarse part of the sponge supplied to remove small impurities from the bottom of the chamber. If lime deposits are found, it is advisable to check the quality of water being used (paragraph 9.5).

QUARTERLY MAINTENANCE 14.3

Quarterly maintenance requires lubrication of the door hinges and replacement of the bacterial filter.

- **HINGE LUBRICATION:** spray the two door hinges with small quantities of silicone oil.
- **REPLACE THE BACTERIAL FILTER:** replace the filter approximately every 300 cycles (NB: bacterial filter life largely depends on usage rather than time. Notwithstanding this fact, it is recommended that the filter should be changed at least every three months, so as not to reduce the autoclave drying performance). A clogged bacterial filter may lead to some difficulties in opening the door on completion of the sterilisation cycle.

YEARLY MAINTENANCE 14.4


The sterilizer is fundamental to the protection of patient and operator alike: even though the electronic controls of these machines are increasingly reliable, it is good practice to carry out a functional check of the equipment at least once a year. This check should be carried out by authorised service centres only, with calibrated and certified instruments, in order to guarantee the equipment a long life and reliable operation (validation). To establish the proper check procedures, refer to the instructions issued by the relevant health authority.

- **YEARLY VALIDATION:** validation requires the use of instruments calibrated by specialist centres to check the sterilisation cycle parameters. This check includes inspections of pressure and temperature probes and the timer. On request, the manufacturer issues a yearly test certificate for machines returned to its premises for maintenance and checks.


14.5 REPLACING THE DOOR GASKET

To replace the door gasket:

- Open the door, use a screwdriver to unscrew the four screws and remove the external plastic/aluminium guard. Press lightly on the dome inside the door in the direction shown by the arrow (fig. A) and at the same time loosen the bolt on the outside of the door (fig. B), without unscrewing it completely.
- Remove the gasket completely from its housing by pulling lightly outwards (fig. C).

- Insert the new gasket. The edge of the steel dome must fit perfectly in the gasket housing (fig. D).
- Press lightly on the dome (fig. E) to prevent it turning and screw up the bolt loosened previously (fig. F). During this operation, make sure the gasket enters the door housing correctly.

- Screw up the bolt until it is completely aligned with the internal screw (fig. G). Block the bolt in position by placing silicon on the external surface.
- Make sure the dome moves slightly along the vertical (fig. H) and horizontal axes (fig. I) by pressing it at either end alternately. This allows it to adapt to the surface of the boiler edge once the door is closed.


Fig. G

Fig. H


Fig. I

- Reposition the external guard and make sure the door closes correctly.

UNSCHEDULED MAINTENANCE 14.6

Unscheduled maintenance (to be carried out by qualified, authorised personnel only) includes repairs of accidental breakages and replacement of worn or faulty parts.

All electrical and mechanical work must only be carried out by an authorised Zhermack maintenance technician.

Modifications which alter the machine characteristics from a safety and risk prevention viewpoint may only be made by the Manufacturer who shall then issue a certificate stating that the machine conforms to safety standards. Therefore any modification or maintenance tasks not actually described in this manual are to be considered forbidden. If the machine is modified, tampered with or repaired by unauthorised personnel or non-original spare parts are used the manufacturer is relieved of all responsibility. If a maintenance task is required but is not described in this manual you must contact an Authorised Servicing Centre.

RESTARTING AFTER A LONG PERIOD OF INACTIVITY 14.7

Should the Zetaclave B need to be restarted after a particularly long period of inactivity proceed as described in Chapter 6 (First-time Use).

If there are any faults or malfunctions proceed as described in Chapter 13 (Alarms and errors). If the problem persists or is not covered by the information in this chapter contact your Authorised Servicing Centre immediately.

REQUESTING TECHNICAL ASSISTANCE 14.8

To obtain technical assistance choose one of the following options:

1. contact your Authorised Servicing Centre;
2. telephone the Zhermack Customer Assistance Service on the number **+39 0425.597.611** and proceed as follows: give personal data (e.g. address and telephone n.), state the information on the equipment ID plate (see paragraph 3.4) and explain the problem clearly;
3. fax all the above information and a properly filled out ID Form (see Appendix A.2) to **+39 0425.53.596**.

APPENDICES

CE DECLARATION OF CONFORMITY

A.1

WARRANTY CERTIFICATE AND IDENTIFICATION FORM

A.2

MAINTENANCE AND REPAIR WORKSHEET

A.3

AUTHORIZED SERVICE PARTNERS

A.4

NOTES

A.5

DECLARATION OF CONFORMITY CE A.1

The Manufacturer:

Person authorised to compile the technical file:

Zhermack S.p.A.
Via Bovazecchino, 100
45021 Badia Polesine • RO • Italy
Tel. +39 0425 597 611 • Fax +39 0425 53 596

Zhermack S.p.A.
Via Bovazecchino, 100
45021 Badia Polesine • RO • Italy
Tel. +39 0425 597 611 • Fax +39 0425 53 596

declares under its sole responsibility that the product:

Model / Serial Number	Lt.
Zetaclave B	18 ☐
	23 ☐

is in conformity with Directive 2011/65 EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

It is also in conformity with the laws implementing the following directives:

- EC Directive 2004/108 EC (EMC Directive) and subsequent amendments.
- EC Directive 2006/95 (Low Voltage Directive) and subsequent amendments.

Medical devices:

EC Directive 93/42 - class IIb

Pressure tool:

EC Directive 97/23.

Zhermack S.p.A. prohibits the use of the medical device referred to in this declaration for purposes other than those indicated in the INSTRUCTION MANUAL.

The device is designed and manufactured in compliance with the basic requirements prescribed in annex I to Legislative Decree no. 46/97 (annex I to Directive 93/42/EEC). When designing and manufacturing this device, best practice standards have been followed:

UNI CEI EN ISO 14971:2012

Medical devices: application of risk management to medical devices.

NOTIFIED BODY

TÜV SÜD (n° 0123). The process chosen for the Evaluation of conformity is as per Annex II to Directive 93/42/EEC.

Please note that the manufacturer's liability for damages caused by a defective product expires 10 years from the date on which the manufacturer released such product on the market and that such liability is transferred to the user in accordance with Directive 85/374/EEC, subsequent amendments and the laws implementing it.

Badia Polesine, 25/02/2014

Paolo Ambrosini
General Manager

A.2 WARRANTY CERTIFICATE AND IDENTIFICATION FORM

1. This document certifies that the Manufacturer has built the product correctly using only top quality materials, carried out all the necessary tests and that the product complies with EC directives. **This product Warranty is valid 24 months or 1500 cycles for the whole product, and 5 years or 5000 cycles for boiler alone, as from date of delivery to the customer. Proof of such date shall be provided by a purchase document or the completed identification form. Labour is excluded from the warranty when not provided by the manufacturer.** These documents must be handed over together with the equipment to the Servicing Centre. Persons wishing to make use of warranty coverage must make the fault known within eight days of its discovery as per Italian civil code art. 1495. Warranty cover is limited to replacement or repair of individual components or pieces with manufacturing defects and does not cover costs related to transfer of technical personnel, transport, packing etc. Breakdowns, faults or damages which derive from poor maintenance, incorrect power supply, negligence, inexperience and other causes not imputable to the Manufacturer and parts subject to normal wear are not covered by the Warranty. The Warranty does not cover breakdowns caused by failure to carry out routine maintenance or customer negligence. This Warranty does not provide for any awards for direct or indirect damages of any kind caused to persons or things by any equipment inefficiency.
2. The Warranty is instantly rendered null and void if the equipment is repaired, modified or tampered with by the purchaser or any unauthorised third parties. Responsibility for damages caused by defective products expires 10 years after the Manufacturer puts them on the market: after such time damages are the user's responsibility, as stated in EC Directive 85/374.
3. To have work carried out under Warranty the purchaser must only contact the retailer, the Manufacturer-authorised servicing centres, or the Manufacturer. The Warranty entitles the purchaser to free replacement of defective parts. It does not entitle the purchaser to replacement of the entire equipment.
4. In the event of disputes over application of the Warranty, the quality or condition of the equipment, the purchaser shall not suspend payment or delay payment of the price or price instalment.
5. The purchaser shall not be entitled to damages for any loss of production.
6. The Warranty is rendered null and void if:
- a. the equipment is damaged by a fall, exposed to fire, spilt liquids, lightning, natural disaster or any other cause not imputable to the Manufacturer.
 - b. installation has not been carried out properly.
 - c. power supply is incorrect (e.g. rated mains voltage incorrect). Erroneous mains connection (rated mains voltage incorrect), or failure to install proper protection devices.
 - d. the serial number has been removed, cancelled out or altered.
7. Any components to be replaced under the terms of the Warranty shall be sent back to Zhermack which shall then send the appropriate spare parts. Should the changed part not be replaced it shall be charged to the person who ordered it.
8. The Manufacturer and the Dental Deposit are not required to provide replacement equipments while repairs are carried out.
9. Spare parts shall, for tax purposes, only be provided under warranty where the relevant limitations set out in the terms and conditions of the warranty are observed.
10. For any other eventuality not provided for by this Warranty Certificate and the regulations consult the Italian Civil Code.
11. Labour, travel and call-out invoices must be paid for on their receipt.
12. The manufacturer and the dealer shall comply with all data processing laws, including proper observance of data security standards in compliance with the specifications given in the data processing information sheet.

IDENTIFICATION FORM

EQUIPMENT MODEL:

EQUIPMENT SERIAL NUMBER:

PURCHASER:

ADDRESS:

TOWN / POSTCODE:

TEL:

TAX CODE / VAT NUMBER RETAILER:

RETAILER:

PURCHASER'S SIGNATURE:

MAINTENANCE AND REPAIR WORKSHEET **A.3**[illegible]

A.4 AUTHORIZED SERVICE PARTNERS

Please visit the page:

http://www.zhermack.com/Technical/Equipment/Authorized_service_partners.kl

to check all the updated telephone numbers and addresses of the Authorized service partners.

A.5 NOTES

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and extend across the width of the page. There is no handwriting or other markings on the paper.

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