**EN - DANGER:** Harmful if swallowed. Causes serious eye damage. Causes skin irritation. May cause an allergic skin reaction. Wash hands thoroughly after handling. Wear protective gloves/eye protection/face protection. IF SWALLOWED: Call a POISON CENTER or if you feel unwell. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Contains: MALEIC ACID, DISODIUM CARBONATE, COMPOUND WITH HYDROGEN PEROXIDE (2-3).

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specific for instruments


concentrated powder | use at 2%

- Bactericidal
- Yeasticidal
- Mycobactericidal, including tuberculocidal
- Virucidal
- Sporidical

Last Update: 2022-03

0118060-2022-03

900 g | sterilizer for medical devices | www.zhermack.com

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**EN** - Cold chemical sterilizer and high-level powder disinfectant for ordinary, surgical and rotating instruments. Zeta 2 Sporex is ideal for all medical devices that can not be sterilized by autoclaving. The product is not suitable for objects made of copper, brass, aluminium alloys, nickel alloys or damaged nickel-chromium plated instruments.

### Use

Zeta 2 Sporex is used as a 2% solution. For every litre of water, add 20 g (3 scoops) of Zeta 2 Sporex. Stir to facilitate dissolution of the powder. Wait for activation of the solution: 15 minutes with water temperature of about 35°C (warm water), at least one hour with water temperature below 35°C (room temperature water). Soak the instruments for 10 minutes. Remove the instruments from the bath, rinse and dry them. We recommend following the recommended soak time. The prepared solution is stable for at least 24 hours but should be renewed at the start of each work day.

Range of action:

Micro organisms	EU standards	Concentration	Reaction time (minutes)
-----------------	--------------	---------------	-------------------------

Bactericidal (S. aureus, P. aeruginosa, E. hirae)	EN 13727 EN 14561	0.5% 0.5%	5 5
Yeasticidal (C. albicans)	EN 13624 EN 14562	0.5% 0.5%	5 5
Mycobactericidal, including tuberculocidal (M. terrae, M. avium)	EN 14348 EN 14563	1% 1%	10 10
Virucidal (tested on Polioviruses, Adenoviruses, Norovirus, Parvovirus including HIV, HBV, HCV)	EN 14476	2%	1
Sporicidal (B. subtilis)	EN 13704	2% 0.5%	5 10

### Composition

Oxygen based (> 40%) bleach; activator (> 30%), anionic and non-ionic surfactants (<5%), corrosion inhibitors, additives and auxiliaries. Dissolution of the powder in water leads to the development of peracetic acid.

### Warnings, storage conditions and precautions

The product is intended for use by qualified operators, with knowledge of the required safety standards. Use protective gloves and appropriate PPE. Operate in well-ventilated areas. Keep the product tightly closed in a dry and cool place. **IMPORTANT:** the undissolved powder remaining on the bottom guarantees the effectiveness of the solution for its entire duration in accordance with the instructions for use. Rinse and clean the dilution pan before introducing Zeta 2 Sporex. Do not mix Zeta 2 Sporex with any other disinfectants.

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**ZHERMAK S.p.A.** Via Bovazecchino, 100, 45021 Badia Polesine (Ro), Italy