

## hydrorise implant Ouick

## minutes

# hydrorise implant

### High-rigidity A-Silicone designed for implantology

#### **OVERVIEW**

- Only 2 minutes in the mouth
- Total setting time of 3 minutes 30 seconds
- Hardness 60 Shore A
- Available in medium consistency for monophase technique
- 380 ml cartridges compatible with the most common 5:1 automatic mixers

Starting from now, you do not have to seek the compromise between speed and precision. The new Quick Set version of Hydrorise Implant achieves the same precision as ever, but faster.

It is the ideal solution for professionals looking for a shorter working time and quicker time in mouth compared to the normal version, thus delivering more efficiency to the dental practice.

Get the same accuracy you are used to in only 2 minutes in the mouth!

#### SAME PRECISION, LESS TIME.



The ideal combination for an **accurate impression** in implantology.

> The **high rigidity** helps holding the transfer in the three-dimensional position.

At the same time, the **high** elastic recovery minimizes the risk of deformation during removal from the mouth. Achieve **more**, in **less** time.

> With a **time in mouth of** only 2 minutes, it speeds up the procedure and make your patient feel more comfortable.

> > minutes

Open to new technologies. Open to innovation.

> Hydrorise Implant, Quick version, enters the digital workflow to improve your communication with the dental lab.

The impression is scannable\* with the most common scanning systems in the market without using matting sprays and also CBCT.

\*Scanning test performed by 3Dfast on impressions obtained from models.



Reliability and safety.

> It's practical! It's safe! Hydrorise Implant is radiopaque, so residues under the gum can be easily detected though X-ray.

**But there's more!** Being an A-Silicone it is **biocompatible** for a greater patient safety.

#### Fulfilling your needs

Zhermack S.p.A. | Via Bovazecchino, 100 45021 Badia Polesine (RO) Italy T +39 0425 597611 | F +39 0425 53596 info@zhermack.com | www.zhermack.com

