



occlufast+
occlufast+ color

occlufast+ occlufast+ color

A-Silicones for bite registration

control
comfort
time saving
precision

Occlufast+ is the new line of A-Silicones for a precise and accurate bite registration.

The imperceptible consistency, that allows to not interfere with the occlusion¹, combined with a time in mouth of only 1 minute, offers a greater comfort to the patient and helps to ensure reliable results.

Your key to a restoration that fits!

Quick **processing**.

MAKE IT FASTER.

Thanks to its **quick time in mouth**, Occlufast+ increases patient comfort and allows the dentist to complete the procedure in a quicker way.

5
minute

A reliable **registration**.

MAKE IT ACCURATE.

Occlufast+ has a **high final hardness** (37 ± 5 Shore D) that helps keeping an accurate and reliable bite registration.

Smart **extrusion**.

MAKE IT SMART.

The easy handling, due to **ready-to-use mixing cartridges**, allows to directly apply the material in the mouth of the patient comfortably. Additionally, the use of the **yellow mixing tip** saves material at each use if compared with other tips*.

Your **visual guide**.

MAKE IT EASIER.

Occlufast+ Color is **thermochromic** and changes color according to the temperature (from green to yellow). The **color change indicates the completed setting** of the material offering a visual guide to the dentist in order to minimize the risk of errors.

¹N. Mobilio & S. Catapano, *Effect of experimental jaw muscle pain on occlusal contacts*,
Journal of Oral Rehabilitation, 2011.

F120052 - 19-05

Fulfilling your needs

Application pictures courtesy of: Dr. A. Barbaglia

Zhermack S.p.A. | Via Bovazecchino, 100
45021 Badia Polesine (RO) Italy
T +39 0425 597611 | F +39 0425 53596
info@zhermack.com | www.zhermack.com

Zhermack 