

cyclope



Equipment
Electronic bunsen unit

Zhermack
Dental

ELECTRONIC INDUCTION TOOL HEATER

cyclope

SAFETY

- No risk of burns for the operator*
- Better air quality in the dental laboratory thanks to elimination of the naked flames which remove oxygen from the work environment

TIME AND MONEY SAVINGS

- The instrument heats up in just 2"
- Energy saving thanks to two external sensors which activate heating only when there is an instrument in the chamber

EASE OF USE

- **Mobile and ready-for-use, cyclope** does not require connection to the main gas supply and can be placed in the most convenient position for the operator
- **Easy-to-read display of the various work phases** with two LEDs
- **Quick cleaning of the chamber** from wax residues thanks to the silicone caps provided

* Scrupulously follow the instructions in the operator and maintenance manual.



Technical data

Code	Model	Dimension	Weight	Mains voltage (V - Hz)	Absorbed power
C308300	Cyclope 230V	L9,1 / P20,9 / H8,9 cm	~ 0,5 kg	230V - ± 10% 50/60	100 Watt

Standard equipment

Code	
XR0100000	Protective caps

For further information on this product, visit the zhermack.com website



BRANCHES

GERMANY
Zhermack GmbH Deutschland
Öhlmühle 10 - D-49448 Marl

T +49 (0) 54 43 / 20 33 - 0
F +49 (0) 54 43 / 20 33 - 11
info@zhermack.de
www.zhermack.de

POLAND
Zhermapol Sp. z o.o.
ul. Augustówka 14
02 - 981 Warszawa - Polska

T +48 22 858 82 72
+48 22 292 71 17
F +48 22 642 07 14
biuro@zhermack.com
www.zhermapol.pl

USA
Zhermack Inc.
P.O. Box 4195, River Edge
New Jersey 07661 - 4195

T +1 (877) 819 6206
+1 (732) 389 8540
F +1 (732) 389 8543
info@zhermackusa.com
www.zhermackusa.com

Zhermack S.p.A. | Via Bovazecchino, 100 | 45021 Badia Polesine (RO) ITALY
T +39 0425 597611 | F +39 0425 53596 | info@zhermack.com | www.zhermack.com

Zhermack