

For further information on these products please visit the website [en.zhermack.com](http://en.zhermack.com)



**zetalabor** **TOUCH**  
**platinum85**



## **FLAVOURED SMOOTH PRECISION**

**Touch the next generation lab putty**


From the inventors of zetalabor, here is the new lab putty **Platinum 85 TOUCH**.

With a smooth and velvety surface, high detail definition and vanilla flavour: discover the characteristics that make this A-Silicone the ideal choice for all applications that require quality and reliability.

**Platinum 85 TOUCH:** feel the next generation Zhermack lab putty.



# A-Silicone lab putty



## Benefits

Smooth surface after polymerization

High detail definition

Dimensional stability over time

Optimised working times

Flavoured

Resistant to high temperatures

## APPLICATIONS

### MASKS for:

- Making temporary prosthetic works
- Repairing complete/framework prostheses
- Creating gingival masks on the model
- Designing metal structures
- Duplication of total removable prostheses and gypsum models
- Verticulator moulding

### BASES for:

- Blocking of prostheses during repairs
- Light curing temporaries in the flask

### MOULDS for:

- Isolation of teeth in the flask

### BITES for:

- Occlusal registration of model in the articulators

### Technical features

Mixing ratio	Flavour	Working time*	Removal time (min.)*	Hardness	Colour	Specific weight
1:1	Vanilla	1'	4'	85 Sh-A	Dark Blue	1.6 g/ml

### Packaging

Code	Description
<b>C400750</b>	1x 4.3 Kg (2690 ml) tub Base + 1x 4.3 Kg (2690 ml) tub Catalyst
<b>C400751</b>	1x 200 g (125 ml) tub Base + 1x 200 g (125 ml) tub Catalyst

\*The specified times may vary depending on the ambient temperature and the mixing technique

### BRANCHES

#### GERMANY

**Zhermack GmbH Deutschland**  
Öhlmühle 10 - D-49448 Marl

T +49 (0) 54 43 / 20 33 - 0  
F +49 (0) 54 43 / 20 33 - 11  
info@zhermack.de  
www.zhermack.de

#### POLAND

**Zhermapol Sp. z o.o.**  
ul. Augustówka 14  
02 - 981 Warszawa - Polska

T +48 22 858 82 72  
+48 22 292 71 17  
F +48 22 642 07 14  
biuro@zhermack.com  
www.zhermapol.pl

#### USA

**Zhermack Inc.**  
P.O. Box 4195, River Edge  
New Jersey 07661 - 4195

T +1 (877) 819 6206  
+1 (732) 389 8540  
F +1 (732) 389 8543  
info@zhermackusa.com  
www.zhermackusa.com