

sympress



As easy as pressing a button

Sympress is an **automatic mixing machine** specifically designed for a **fast** and **easy mixing** of **impression materials packaged in 5:1 cartridges**.

Mixing becomes automatic to ease your daily work.

VERSATILE, RELIABLE, EASY TO USE.

Zhermack 
Dental



Versatile
Open system



Reliable
Rely on quality



Easy to use
Work smarter

- ▶ Compatible with the most common impression materials in the market (VPS and polyether)
- ▶ Can be used both with hard and soft cartridges (360 or 380 ml)
- ▶ German quality you can trust*
- ▶ Functionality: compact design for reduced impact on counter space
- ▶ Low noise
- ▶ Easy to use to offer maximum flexibility and comfort to dental clinics
- ▶ Quicker processing to increase efficiency
- ▶ Mixing is uniform and delivers a higher quality than manual mixing

Try Sympress with our 5:1 systems

FILL LEVEL INDICATOR

showing the amount of material remaining in the cartridge

LED CONTROL PANEL

indicating if the device is ready

TWO DISPENSING SPEEDS

fast for tray filling and slower for syringe filling

AUTOMATED LOADING

for easy setup

SHAPED TRANSPARENT COVER

for the correct fitting of different types of cartridges

SMALL FOOTPRINT

with the flexibility of wall mounting



380 ml cartridges for an homogeneous and efficient mixing.



Patented tip

Reduced material waste of 18% in comparison to the traditional mixing tips.**

* Sympress is offered in partnership with Renfert. ** Zhermack internal tests.

Technical features

Code	Model	Size (H x W x L - mm)	Weight (kg)	Power supply (V/Hz)	Total power (W)
6000-0000	Sympress - 230 V	330 x 190 x 380	7	230 V ± 10% / 50 Hz	170
6000-1000	Sympress - 120 V	330 x 190 x 380	7	120 V ± 10% / 60 Hz	170
6000-2000	Sympress 230 V plug UK	330 x 190 x 380	7	230 V ± 10% / 50 Hz	170
6000-3000	Sympress 100 V	330 x 190 x 380	7	100 V ± 10% / 50 - 60 Hz	170

Accessories

Code	Description
6000-0100	Wall mounting kit