



zeta 1 ultra

Zeta 1 Ultra is a concentrated liquid disinfectant and detergent with a new, efficient formula for ordinary and surgical dental instruments and burs.

This new and highly effective product can be used in two ways:

- **Broad spectrum of action:**
use at **2%** for 60 minutes in a static bath, or for 30 minutes at 35°C in an ultrasonic bath for complete virucidal action.
- **Limited spectrum of action:**
use at **1%** for 15 minutes in a static bath for a bactericidal, yeasticidal and limited virucidal action (HIV, HBV, HCV).

Zeta 1 Ultra, safe and effective.

www.zhermack.com



Zhermack 
Dental